# tonocare

## Tonometer

INSTRUCTIONS FOR USE

EN DE ES FR IT PT

 Keeler




EN	TonoCare Instructions For Use .....	1
DE	TonoCare Gebrauchsanleitung.....	32
ES	TonoCare Instrucciones De Uso.....	63
FR	TonoCare Mode D'emploi .....	94
IT	TonoCare Istruzioni Per L'uso .....	125
PT	TonoCare Instruções De Utilização .....	156

## CONTENTS

<b>1. INDICATIONS FOR USE</b> .....	<b>5</b>
1.1 INTENDED USE / PRINCIPLE OF OPERATION.....	5
<b>2. SAFETY</b> .....	<b>5</b>
2.1 PHOTOTOXICITY .....	5
2.2 WARNINGS AND CAUTIONS.....	6
2.3 CONTRAINDICATION.....	8
<b>3. CLEANING INSTRUCTIONS</b> .....	<b>8</b>
3.1 CLEAN THE PUFF TUBE LENS WEEKLY: .....	8
3.2 CLEAN THE HAND UNIT DAILY AND BETWEEN PATIENTS .....	9
<b>4. INSTALLATION</b> .....	<b>9</b>
4.1 TOOLS REQUIRED.....	9
4.2 PACKAGING CONTENT .....	9
4.3 TABLETOP INSTALLATION .....	10
4.4 WALL MOUNT INSTALLATION .....	10
<b>5. USING THE TONOCARE</b> .....	<b>11</b>
5.1 CONTROLS AND INDICATORS.....	11
5.2 MENU MAP.....	13
<b>6. MEASUREMENT PROCEDURE</b> .....	<b>14</b>
6.1 ROUTINE CHECKS AND FUNCTIONAL TESTS.....	14
6.2 PREPARING THE DEVICE.....	15
6.3 PREPARING THE PATIENT.....	15
6.4 TAKING THE READING .....	15
6.5 CCT CORRECTION.....	17
<b>7. PRINTING</b> .....	<b>19</b>
<b>8. REPLACING THE PRINTER PAPER</b> .....	<b>19</b>
8.1 CHARGING YOUR TONOCARE.....	20
<b>9. MAINTENANCE</b> .....	<b>20</b>
9.1 FAULT CODES .....	21
<b>10. WARRANTY</b> .....	<b>21</b>
<b>11. SPECIFICATIONS AND ELECTRICAL RATINGS</b> .....	<b>22</b>
11.1 ELECTROMAGNETIC EMISSIONS .....	22
11.2 ELECTROMAGNETIC IMMUNITY .....	22
11.3 RECOMMENDED SAFE DISTANCES .....	24
<b>12. TECHNICAL SPECIFICATIONS</b> .....	<b>25</b>
<b>13. ACCESSORIES AND SPARES</b> .....	<b>27</b>
<b>14. PACKAGING AND DISPOSAL INFORMATION</b> .....	<b>27</b>

<b>15. APPENDIX.....</b>	<b>28</b>
15.1 REPEATABILITY AND REPRODUCIBILITY.....	28
15.2 CLINICAL PERFORMANCE DATA.....	28

	Consult instructions for use		General warning sign
	Date of manufacture		Warning: Non-ionizing radiation
	Manufacturer's name and address		Warning: Electricity
	Country of manufacture		Warning: Optical radiation
	Waste Electrical and Electronic Equipment (WEEE) recycling		Warning: Floor level obstacle
	This way up		Keep dry
	Type BF applied part		Fragile
	Temperature limit		Do not use if package is damaged
	Authorised representative in the European Community		Authorised representative in Switzerland
	United Kingdom Conformity Assessed, with the Notified Body number for SGS UK		Conformité Européenne, with the Notified Body number for SGS Belgium NV
	Catalogue number		Class II equipment
	Serial number		Atmospheric pressure limitation
	Medical device		Humidity limitation
	Translation		

The Keeler TonoCare is designed and built-in conformity with Directive 93/42/EEC, Regulation (EU) 2017/745 and ISO 13485 Medical Devices Quality Management Systems.

Classification: CE / UKCA: Class IIa  
FDA: Class II

The information contained within this manual must not be reproduced in whole or part without the manufacturer's prior written approval. As part of our policy for continued product development we the manufacturer reserve the right to make changes to specifications and other information contained in this document without prior notice.

This IFU is also available on the Keeler UK and Keeler USA websites.

Copyright © Keeler Limited 2023. Published in the UK 2023.

## Symbols used on controls and display

These symbols will be used on the device and display during operation.

Device Controls		Display Symbols		Description
	On / Standby button		Light	Turn On / Off illumination LEDs
	Menu button		Sound	Turn On / Off audible alerts
	Return button		Time	Set the system time
	Manual trigger		Date	Set the system date
	OD / OS button		IOP Format	Select the IOP format (XX / XX.X)
	Print button		Pachymetry	Pachymetry options of OFF/ON/ TRIGGER  (If set to 'TRIGGER', the pachymetry option only appears if the measured IOP value is greater than 15)
<b>Sounds</b>				
<b>LOW PITCH</b>	Bad or invalid readings		Brightness	Set the brightness for the displays (value between 1 and 15)
<b>HIGH PITCH</b>	Valid set of measurements obtained			
			Self-Test	Initiate a self-test of the system

## 1. INDICATIONS FOR USE

These devices are intended to be used only by suitably trained and authorised healthcare professionals.

**The TonoCare Non-contact Tonometer should be used only by trained personnel. USA Federal law restricts this device to sale by or order of a physician.**

### 1.1 INTENDED USE / PRINCIPLE OF OPERATION

The Keeler TonoCare Tonometer is a hand-held, battery operated, non-contact tonometer intended to be used for measuring intraocular pressure (IOP) of the human eye with less than 3D in corneal astigmatism.

The Keeler TonoCare should not be used in patients with high corneal astigmatism (>3D).

The Keeler TonoCare Tonometer uses the principle of air impulse tonometry – this is a variation of general applanation tonometry, however, does not require direct contact with the surface of the eye.

Applanation tonometry is a technology to accurately measure IOP as an equivalent of the force required to flatten a defined area of the cornea by mechanical stimuli, as a direct application of the Imbert-Fick law.

The air impulse technique requires direction of a packet of air with restricted pressure and volume towards the central portion of the cornea, and the detection of the pre-defined flattening of the cornea via the electrical measurement of a light beam reflected from the corneal surface.

## 2. SAFETY

### 2.1 PHOTOTOXICITY

**CAUTION: The light emitted from this instrument is potentially hazardous. The longer the duration of exposure, the greater the risk of ocular damage.**

**While no acute optical radiation hazards have been identified for Keeler Tonometers, we recommend keeping the intensity of the light reaching the patient's retina to the minimum possible for the respective diagnosis. Children, people with aphakia and people suffering from eye conditions are most at risk. An increased risk may also occur if the retina is exposed to the same or a similar device with a visible light source within 24 hours. This applies, in particular, if the retina has been photographed with a flashbulb in advance.**

**Keeler Ltd shall on request, provide the user with a graph showing the relative spectral output of the instrument.**

## 2.2 WARNINGS AND CAUTIONS

Please note that the proper and safe functioning of our instruments is only guaranteed if both the instruments and their accessories are exclusively from Keeler Ltd. The use of other accessories may result in increased electromagnetic emissions or reduced electromagnetic immunity of the device and may lead to incorrect operation.

Observe the following precautions in order to ensure safe operation of the instruments.

### WARNINGS

- Never use the instrument if visibly damaged and periodically inspect it for signs of damage or misuse.
- Check your Keeler product for signs of transport / storage damage prior to use.
- US Federal Law restricts this device to sale by or on the order of a physician or practitioner.
- The device is intended for use in various clinical settings such as hospitals, eye clinics and optometric practices.
- Only use approved Keeler power supply EP29-32777 or instrument may malfunction.
- If transporting the TonoCare in temperatures of less than 10°C, ensure the device has had at least three hours to acclimatise to room temperature before use.
- The owner of the instrument is responsible for training personnel in its correct use.
- Never use the instrument if the ambient temperature, atmospheric pressure, and / or relative humidity are outside the limits specified in this manual.
- Do not use in the presence of flammable gases / liquids, or in an oxygen rich environment.
- This device is intended to be used only by suitably trained and authorised healthcare professionals.
- This product should not be immersed in fluid.
- Do not mount anything on the Docking Station other than the Hand Unit and Forehead Rest.
- TonoCare is not intended to be used with wireless technology. Do not plug a wireless dongle into the USB port on the Docking Station.
- The mains plug is the means of isolating the device from the mains supply. Ensure both the power switch and mains plug are always accessible.
- Do not position the equipment so that is difficult to remove the mains plug from the wall socket.

- Do not fit mains power adapter into a damaged mains outlet socket.

- Using the device out of specification or where there is a fault can lead to an electric shock, overheating and/or burns.

- Route power cords safely to eliminate risk of tripping or damage to user.

**CAUTION**

- Use only genuine Keeler approved parts and accessories or device safety and performance may be compromised.
- Keep out of the reach of children.
- Accuracy of IOP measurements is known to be affected by variations and changes in corneal rigidity due to differences in corneal thickness, intrinsic structural factors or corneal refractive surgery. It is recommended that these factors are considered during IOP measurement.
- To prevent condensation from forming, allow instrument to come to room temperature before use.
- Please add mounting hole labels to covering the key holes unless wall mounting the TonoCare.
- Only mount on wall according to Keeler Instructions.
- The TonoCare cannot be used in the vicinity of sources known to cause electromagnetic disturbance (magnetic resonance imaging, computed tomography, radio-frequency identification, metal detectors, electronic article surveillance and other electromagnetic security systems).
- Keep the front window and nozzle area away from large amounts of dust or fine particles.
- This product should be used in a room with low / dimmed lighting.
- The device will require a minimum 12-hour charge before first use.
- The device will require several 12-hour charge cycles before the battery functions optimally.
- Before using the TonoCare, press the Manual fire button to dispel any minute particles of dust or moisture which may have settled whilst the instrument was not in use.
- Check the function of the device in accordance with the instructions in section 6 prior to use on the patient.
- The patient should not be in the proximity of the Docking Station.
- Do not touch the electrical contacts on the Docking Station and the patient simultaneously.
- The Forehead Rest is composed of aluminium, and it is the only part which may touch the patient. The housing of the instrument is made of PC-ABS. Do not touch these parts if you have a known allergy to any of the materials.
- Do not use the Forehead Rest in the Docking Station for any other instrument, or the performance of the product may be compromised.
- Contact of the front window / nozzle area of the TonoCare with the eye of the patient should be avoided. If accidental contact occurs, clean the front window and surrounding area according to the cleaning instructions in section 3.1.
- Always verify that the printout matches the readings on the Hand Unit.
- For indoor use only (protect from moisture).
- There are no user serviceable parts inside. Contact authorised service representative for further information.

- Follow guidance on cleaning / routine maintenance to prevent personal injury / damage to equipment.
- Failure to carry out recommended routine maintenance as per the instructions in this IFU may reduce the operational lifetime of the product.
- If the device is not used regularly it must be recharged for at least 12 hours monthly to ensure optimum battery life.
- At product end of life dispose of in accordance with local environmental guidelines (WEEE).

### Maintenance

- To maintain device performance, and ensure its safety and effectiveness it should be serviced in accordance with the instructions in section 9.
- Only decontaminate / clean in accordance with the instructions given in section 3.1.
- If the device is not used regularly it must be recharged for at least 12 hours monthly to ensure optimum battery life.

### 2.3 CONTRAINDICATION

There is no restriction to patient population this device can be used with other than those outlined in the contraindications stated below.

Accuracy of IOP measurements is known to be affected by variations and changes in corneal rigidity due to differences in corneal thickness, intrinsic structural factors or corneal refractive surgery. It is recommended that these factors are considered during IOP measurement. The Keeler TonoCare should not be used in patients with high corneal astigmatism (>3D).

## 3. CLEANING INSTRUCTIONS

**CAUTION: Only manual non-immersion cleaning as described should be used for this tonometer. Do not autoclave or immerse in cleaning fluids. Always disconnect power supply from source before cleaning.**

### 3.1 CLEAN THE PUFF TUBE LENS WEEKLY:

1. Moisten a cotton bud with de-ionised water / detergent solution (2% detergent by volume).
2. Move the tip of the bud around the lens in a circular motion.
3. After one circle the bud should be discarded to avoid smearing on the lens.
4. Look at the Puff Tube lens from the patient's side, if traces of tear film can still be seen, repeat above steps until clear.

**Note:** Care should be taken not to damage the Puff Tube assembly during cleaning.

**CAUTION: Care should be taken not to damage the Puff Tube assembly during cleaning. Never use a dry cotton bud or tissue to clean the Puff Tube lens. Never use a silicone impregnated cloth or tissue to clean the puff tube lens.**

### 3.2 CLEAN THE HAND UNIT DAILY AND BETWEEN PATIENTS

1. Wipe the external surface with a clean absorbent, non-shedding cloth dampened with de-ionised water / detergent solution (2% detergent by volume) or water / isopropyl alcohol solution (70% IPA by volume). Avoid using de-ionised water / isopropyl alcohol solution with optical surfaces such as the front window. These should only be cleaned with a de-ionised water / detergent solution
2. Ensure that excess solution does not enter the instrument. Use caution to ensure cloth is not saturated with solution.
3. Surfaces must be carefully hand-dried using a clean non-shedding cloth.
4. Safely dispose of used cleaning materials.

## 4. INSTALLATION

This section will instruct how to unpack and prepare the TonoCare for use.

**CAUTION:** When you open the package, check for any external damage or flaws, particularly damage to the case. If you suspect there is something wrong with the tonometer, contact the manufacturer or distributor.

Replace the blanking plate with the appropriate mains plug adapter if required, or use IEC 60320 TYPE 7 connector (not supplied).

### 4.1 TOOLS REQUIRED

- Safety knife

#### For wall mounting:

- Pencil
- Spirit Level
- Electric drill
- PH1 screwdriver

### 4.2 PACKAGING CONTENT

Your TonoCare has been supplied with:

- A **Hand Unit** with pre-installed battery for IOP standalone measurement.
- A **Docking Station** for printing, data export and charging functions.
- A **metal wall mounting plate** with 4 screws and 4 rawlplugs for wall mounting of the Docking Station.
- A **power supply** for directly charging the Hand Unit (in transport) or via the Docking Station.
- A **roll of thermal paper** to be used in the printer located in the Docking Station.
- A **USB device** containing the Instructions for Use.

- An extendable **Forehead Rest**.
- A **USB cable** to connect the Docking Station with a computer (not provided) for data export.
- A **Tonometer Face Shield** to provide a layer of protection between the patient and user.

#### 4.3 TABLETOP INSTALLATION

1. Position the product package so that the arrow points upward
2. Use a safety knife to cut the tape sealing the box on the top and remove the polystyrene layer covering the package contents.

**Use caution to avoid injury from sharp edges when handling the safety knife, and the unsealed carton edges.**

3. Lift the Docking Station from the package and place it in a clean area intended for the charging of the TonoCare when not in use.

**The TonoCare Hand Unit should not be used in the proximity of the Docking Station.**

4. Insert the roll of thermal paper provided, according to instructions on section 8.
5. Remove the Hand Unit from the package and allow the handle to fit the lower recess in the Docking Station and the Measurement Window to slide onto the top of the Docking Station. Before use, remove the protective film covering the TonoCare display, the Front Window and the IR Window.
6. Take the Forehead Rest from the package and allow it to be held by its magnets at the top of the Docking Station, in the area provided for it.
7. Remove the power supply from the package, plug it to the back of the Docking Station and after fitting the suitable adapter for your country, connect it to the AC inlet.

**Route power cords safely to eliminate risk of tripping to user, or damage to device.**

8. The LED indicator in the Hand Unit should now light up to indicate the TonoCare battery is charging.
9. Use the USB cable provided to connect the Docking Station to a computer (not provided) for data export. The computer must be compliant to EN 60601-1 (see section 12)

#### 4.4 WALL MOUNT INSTALLATION

1. Choose carefully the intended location for your TonoCare Docking Station with particular consideration to the routing of the power cable and the patient screening position. Ensure that the mains plug is accessible at all times, as this is the primary means of mains power disconnection.

The TonoCare Hand Unit should not be used in the proximity of the Docking Station. The Docking Station should not be fixed over live utilities, as the drilling process could interrupt the utility supply and cause injury. The recommended height is 1.2m (4 feet).

- Use the metal plate as a template to mark the position of the retaining screws with a pencil, holding a spirit level to the base of the plate to guarantee horizontal alignment.
- Drill the appropriate size holes following the marks left from the previous step.

Use extreme caution when operating the drill, following the instructions provided with the instrument.

- Insert the rawlplugs in the holes drilled in the previous step, and fix the metal plate on the wall with the screws provided, using a PH1 screwdriver.
- Position the Docking Station on the wall mounting plate so the 2 retaining pins on the metal plate slide into the holes in the rear of the housing, and the plate supports the unit from below.
- You may now plug the power cable into the Docking Station, and connect it to mains power. The Docking Station will blink twice at power up. Once powered, rest the hand

## 5. USING THE TONOCARE

This section will instruct the user on how to interpret the controls and indicators of the TonoCare, and how to perform the IOP measurement on a patient using the device.


Familiarise yourself with the instructions for the measurement of the IOP using the TonoCare, before using the instrument on a patient.

### 5.1 CONTROLS AND INDICATORS

#### Hand Unit

##### Top view

- Forehead Rest mounting
- Manual Trigger
- OD / OS button
- LED indicator
- Digital Display
- On / Off
- Print button
- Menu button
- Return button

**Note:** The LED indicator on the TonoCare pulses when charging and remains solid when fully charged


### Isometric view from the user and left side

- 10 Forehead Rest
- 11 Eyepiece

### View from underneath the tonometer

- 12 Charging contacts
- 13 Power input


### Patient view

- 14 Puff Tube and window
- 15 LED indicator
- 16 Infrared transmitter

### Docking Station

#### Rear view

- 17 Charging contacts
- 18 Infrared receiver
- 19 Printer door
- 20 Wall mounting holes


**Underside view**

- 21 Power input
- 22 USB interface – Allows tonometer to be connected to a PC for raw data upload purposes.

View to show position  
of Hand Unit with  
Docking Station

**5.2 MENU MAP**

Press the Menu button to open the Software Menu. Once in Software Menu, use the Menu button to change your selection, the Print button to confirm and the Return button to move to the previous menu. Follow the Menu Map in the figure on page 14 for guidance.

Confirm (Print button)


Change selection (Menu button)

Return to previous menu (Return button)

The software revision is displayed during start-up

## 6. MEASUREMENT PROCEDURE

### 6.1 ROUTINE CHECKS AND FUNCTIONAL TESTS

- Before daily use, visually examine the Hand Unit and Docking Station, looking for any signs of obvious damage. If you suspect the tonometer has been subjected to any impacts or moisture ingress do not use the tonometer. Contact Keeler or your local service centre for advice.
- Keeler recommends performing a weekly self-test on the Hand Unit to confirm functionality. This verifies operation of the puff generation system and the pressure sensor. To access the self-test facility, follow these steps:
  - Power the tonometer up, ensuring it is plugged into the mains power supply.
  - Rest the tonometer on its side on a flat surface, such as on a desk. Do not hold the tonometer.
  - Press and hold the Menu button and scroll through until the  icon appears on the display.
  - Using the Print or Return button set this to ON.
  - Press and hold the Menu button again to exit the menu (approximately 2 seconds).
  - The tonometer will puff several times and display a series of numbers as it progresses through the test.
  - Once the self-test is complete the tonometer will indicate that it has passed or failed on the display.

**Remember to prepare the patient before initiating the measurement. A patient's anxiety may delay the measurement and adversely affect its accuracy.**

**A single reading can be misleading as the IOP will vary in response to pulse, respiratory and diurnal fluctuations. Other factors may affect IOP, such as blinking, eye squeezing, fluid intake, physical activity, body position, etc. Up to 4 readings may be required to reduce the impact of these variants to a constant IOP. Keeler recommends using the average of four readings rather than any individual reading.**

The TonoCare software will recognise the readings and give an audible notification when two consecutive readings are within 1mmHg of each other, indicating that further measurements may not be required.

## 6.2 PREPARING THE DEVICE

1. Before lifting the Hand Unit from the Docking Station, check the LED indicator is steady on to ensure full battery capacity. A full battery will last up to 2 days of intensive use. Performance will deteriorate over time.
2. Remove the Hand Unit from the Docking Station and press the Power On button. The Hand Unit will enter Standby Mode if not used for more than 90 seconds.

## 6.3 PREPARING THE PATIENT

Before using the TonoCare tonometer you should make your patient feels at ease and ensure they are in an optimum reading location, preferably with their head supported. This is because apprehension and nervousness may adversely affect the readings obtained. Follow the points outlined below to achieve this:

1. Ensure that the patient is comfortable and in a relaxed position.
2. Ask the patient to remove their contact lenses or spectacles if worn and to breathe normally. The patient's eyes should be fully open and blinking normally throughout the complete measurement procedure.
3. To reassure the patient, you can demonstrate the measurement procedure using the Manual trigger button towards one of the patient's fingers. Re-set (long press OD / OS) the tonometer after the demonstration.

### Before taking a reading, you should:

1. Ask the patient to blink to ensure a good and reflective tear film.
2. Ensure the patient and tonometer optics are not positioned under direct lighting (i.e. spot lights or sunlight).
3. Ensure the patient's eyes are fully opened. This helps to prevent squeezing, where the patient unconsciously tenses their eyelids and increases IOP.
4. Throughout the reading process, you should allow the patient to blink at intervals to maintain the corneal tear film.

## 6.4 TAKING THE READING

Once the TonoCare and the patient are prepared, you are ready to take a reading.

1. If you haven't done it in the previous step, use the manual trigger / demo button to dispel any minute particles of dust or moisture which may have settled whilst the TonoCare was not in use.

**The TonoCare is set to automatically select the right eye as the first eye to be measured. If you wish to select the left eye, press the OD / OS button. Ensure you are recording the measurement for the intended eye.**

2. Hold the Hand Unit with the dominant hand and position the tonometer so it is aligned with the patient's eye from a distance of about 30cm or 12 inches.
3. Move the tonometer in towards the patient until the Forehead Rest is positioned against their forehead. You may choose to rest your fingers on the Forehead Rest to gain more stability.
4. Ask the patient to focus on the green target inside the tonometer.
5. Using your preferred eye, look through the eyepiece so that you can see the patient's eye. Keep moving the tonometer forward whilst aligning the positioning ring to the patient's limbus. The entire iris does not need to be visible as it may be covered by an eyelid. Centre the measurement area with the pupil of the patient's eye and make note of the reflected LED crescents. These should be central on the pupil as well.
6. When the tonometer is approximately 15mm (0.5 inches) away from the patient's eye, a cross will appear to indicate the position of the tonometer relative to the eye. Move the tonometer so that the edges of the cross are just inside the corners of the measurement brackets.

**The cross must not be larger than the measurement area, as this would indicate that the tonometer is too close to the eye.**


7. Once the cross reaches the correct alignment position, a gentle air impulse will trigger an IOP reading. Ensure eyelids and eyelashes are clear of the measurement brackets to obtain accurate results.

**If no applanation was recorded during the puff, a low-pitch tone will be heard (if sounds are enabled in the menu setting) and two stars (\*\*) will be shown on the internal display.**

8. Ensure an IOP reading has been recorded on the tonometer.

**\*Correct position and size of cross in measurement area**


9. Slowly move the tonometer backwards and allow the patient's eye to rest for a few seconds, maintaining the alignment position.
10. When the patient is ready for another reading, move the tonometer closer until the alignment cross appears again and another reading is triggered.
11. Repeat the previous steps for further readings until the measurement averaged from the individual readings is acceptable.
12. When two consecutive readings are within 1mmHg, an audible high pitch tone will be heard indicating that sufficient readings may have been taken (if Buzzer is enabled in the User Menu Options). If successive readings are not within 1mmHg of each other, Keeler recommends taking up to four readings and using the average.

13. Press the OD / OS button to switch from left to right eye, or vice-versa.
14. Press and hold the OD / OS button to clear all the readings.
15. Pressing the Print button will produce a paper printout from the printer installed in the Docking Station. The infra-red window in the Docking Station and Hand Unit should be unobstructed and aligned within 1m (3 feet). The Docking Station LED will flicker during the transmission of the IR data and then extinguish when printing.
16. If the Docking Station is connected to a computer, pressing the Print button will export the raw data to the computer, as long as the serial port has been enabled as described below.

## 6.5 CCT CORRECTION

1. Press and hold the Menu button to open the Software menu. Follow the Menu Map in section 5.2 on page 14 to ensure the CCT correction function is enabled.
2. Follow the instructions in section 6.4 on page 15 to retrieve an IOP measurement. Once a suitable average reading is obtained, follow the CCT workflow shown below. The CCT group selected for the patient should correspond to the one measured separately by a pachymeter. Use the Menu button to confirm your selection.
3. The CCT screen will show the CCT group and the correction applied under the relevant eye.

## TonoCare CCT Workflow

The accuracy of IOP measurements is known to be affected by variations and changes in corneal rigidity due to differences in corneal thickness, intrinsic structural factors or corneal refractive surgery. It is recommended that these factors are considered during IOP measurement. The biomechanical properties of an individual cornea may vary, resulting in changes of the relative stiffness or rigidity of the cornea and altering the measurement. Other factors to consider include corneal edema and other corneal abnormalities potentially affecting rigidity (e.g., keratoconus, corneal transplant, crosslinking) in addition to intrinsic structural factors and corneal refractive surgery.

## 7. PRINTING

The results can be printed by pressing the Print button on the Hand Unit.

This will automatically include the date and time (if set).

A space is included for manually recording the name of the patient.

The last four individual readings are printed as whole numbers 'XX'.

The average IOP is calculated and printed to one decimal place 'XX.X'.

**Always verify that the data on the printout and data on the TonoCare Export application matches the readings on the Hand Unit.**

**The printer on the docking station contains a sharp serrated blade to help cut the paper. Exercise caution to avoid contact with this blade whenever replacing printer paper rolls or tearing printouts from the device.**

The TonoCare Docking Station can also be connected to TonoCare Export through the USB port on a PC. Measurement data can then be transferred from the Hand Unit to the application via the Docking Station.

Refer to EP59-47228 TonoCare Export App Installation Guide for details on how to install the application on your PC. EP59-47333 TonoCare Export App Quick Reference Guide describes all the features of the application and how to use it.

Pressing the Print button will send the result to both the printer and the USB port.

KEELER

Name:  
Date: DD/MM/YY  
Time: HH:MM

Eye:           L       R

Data:           0       0  
                  0       0  
                  0       0  
                  0       0

Avg\_IOP:   0.0   0.0

**Print Example**

## 8. REPLACING THE PRINTER PAPER

1. Access to the printer paper is via the Printer Cover, pull the lip on the top of the Cover and gently pull towards you to open the Printer housing.
2. Remove the empty paper roll.

3. Place the new roll of paper into the paper holder, making sure the free end is loose at the top of the roll and oriented as shown.
4. Extend a few centimetres of paper out of the housing. While holding the end of the paper, close the Cover by gently pushing the lip at the top towards the Docking Station until it is fully closed and clicks into place.

**The printer on the Docking Station contains a sharp serrated blade to help cut the paper. Exercise caution to avoid contact with this blade whenever replacing printer paper rolls or tearing printouts from the device.**

### 8.1 CHARGING YOUR TONOCARE

When not in use Keeler recommends storing your TonoCare on the Docking Station, so that it is maintained fully charged and ready for use.

**The LED on the TonoCare pulses when charging.**

**Once fully charged the LED will be continuously lit.**

The LED on the Docking Station will not change when the TonoCare handset is stored on the Docking Station.

## 9. MAINTENANCE

**Keeler recommends routine maintenance be carried out by the user frequently and an annual service of your TonoCare by Keeler Service Engineers to ensure safe and accurate measurement. In the event of the device being outside of the calibration tolerances, it is important to send the device back to Keeler Ltd. or your local dealer for repair and re-calibration.**

There are no user replaceable parts inside the device, including the battery. The battery must only be replaced by trained service personnel following the instructions in the Service Manual.

If you notice a significant reduction of the battery's performance, contact Keeler or your authorised distributor for its replacement.

If the TonoCare is dropped by accident, the same service centre or distributor can verify if the device is still in calibration.

The unit performs a self-function check when switched on and will indicate if a fault is found. A further self-test can be activated from the menu (refer to section 6.1).

Always inspect the product before use, checking for normal start-up.

Do not attempt to perform product disassembly, reassembly, or repair. These should only be done by personnel trained and qualified by Keeler, following the instructions in the service manual.

Do not store the product in a dusty environment, as the dust may enter the puff system, and be dispelled to the patient's eye during use.

If the TonoCare is to remain unused for any length of time, press the On / Off Push button switch to 'Off' and remove the power supply. Use the dust cover to protect the tonometer.

**Do not attempt to perform any unauthorized repairs, as this could endanger the product and patients. Do not allow unauthorized parts to be fitted in your product.**

On request, Keeler will provide necessary circuit diagrams, component parts lists, descriptions and calibration instructions to assist service personnel in device repair.

The MOD RECORD label on the rear of the device is used to indicate the status of the device in relation to significant changes.

MOD RECORD			
1	6	11	16
2	7	12	17
3	8	13	18
4	9	14	19
5	10	15	20

## 9.1 FAULT CODES

If an Error Code between 00 to 34 is seen on the screen, restart the tonometer and check its function. If the tonometer does not clear, please return it to your nearest authorised Keeler Service Centre.

## 10. WARRANTY

Your Keeler product is guaranteed for 2 years and will be replaced, or repaired free of charge subject to the following:

- Any fault due to faulty manufacture.
- The instrument and accessories have been used in compliance with these instructions.
- Proof of purchase accompanies any claim.

**The manufacturer declines any and all responsibility and warranty coverage should the instrument be tampered with in any manner or should routine maintenance be omitted or performed in manners not in accordance with these manufacturer's instructions.**

**There are no user serviceable parts in this instrument. Any servicing or repairs should only be carried out by Keeler Ltd. or by suitably trained and authorised distributors. Service manuals will be available to authorised Keeler service centres and Keeler trained service personnel.**

## 11. SPECIFICATIONS AND ELECTRICAL RATINGS

The Keeler TonoCare is a medical electrical instrument. The instrument requires special care concerning electromagnetic compatibility (EMC). This Section describes its suitability in terms of electromagnetic compatibility of this instrument. When installing or using this instrument, please read carefully and observe what is described here.

Portable or mobile-type radio frequency communication units may have an adverse effect on this instruments, resulting in malfunctioning.

### 11.1 ELECTROMAGNETIC EMISSIONS

#### Guidance and manufacturer's declaration – electromagnetic emissions

The Keeler TonoCare is intended for use in the electromagnetic environment specified below. The customer or user should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Keeler TonoCare uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Keeler TonoCare is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	N/A, Class A < 75W	
Voltage fluctuations / flicker emissions IEC 61000-3-3	N/A, Class A < 75W	

### 11.2 ELECTROMAGNETIC IMMUNITY

#### Guidance and manufacturer's declaration – electromagnetic immunity

The Keeler TonoCare is intended for use in the electromagnetic environment specified below. The customer or the user should assure that it is used in such an environment.

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD). IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst. IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines 100kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment – guidance
Surge. IEC 61000-4-5	$\pm 1$ kV line(s) to neutral	$\pm 1$ kV line(s) to neutral	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines. IEC 61000-4-11	$<5\% U_T$ ( $>95\%$ dip in $U_T$ ) for 0.5 cycle $<5\% U_T$ ( $>95\%$ dip in $U_T$ ) for 1 cycle $40\% U_T$ (60% dip in $U_T$ ) for 5 cycles $70\% U_T$ (30% dip in $U_T$ ) for 500ms $<5\% U_T$ ( $>95\%$ dip in $U_T$ ) for 5s	$U_T = 0\%$ 0.5 cycle (0, 45, 90, 135, 180, 225, 270, 315°) $U_T = 0\%$ ; 1 cycle and 5 cycles $U_T = 70\%$ ; 25/30 cycles (Single phase: at 0°) $U_T = 0\%$ ; 250/300 cycle	Mains power quality should be that of a typical commercial or hospital environment.  If the user of the Keeler requires continued operation during power mains interruptions, it is recommended that the charger be powered from an uninterruptible power supply.
Power frequency (50/60 Hz) Magnetic field. IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at a level characteristic of a typical location in a typical professional healthcare facility environment.

Note:  $U_T$  is the a. c. mains voltage prior to application of the test level.

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment – guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the Keeler TonoCare, including cables, than the recommended separation distances calculated from the equation applicable to the frequency of the transmitter.
<b>Recommended separation distance</b>			
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 V	$d = 1.2 \sqrt{p}$
Radiated RF IEC 61000-4-3	3 V/m 80MHz to 2.7GHz	3 V/m	$d = 1.2 \sqrt{p}$ 80MHz to 800 MHz $d = 2.3 \sqrt{p}$ 800MHz to 2.7GHz

			<p>Where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey<sup>1</sup>, should be less than the compliance level in each frequency range.<sup>2</sup></p> <p> Interference may occur in the vicinity of equipment marked with this symbol.</p>
--	--	--	---

Note 1: At 80MHz and 800MHz, the higher frequency range applies.

Note 2: These guide lines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people

<sup>1</sup> Field strengths from fixed transmitters, such as base stations (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Keeler TonoCare is used exceeds the applicable RF compliance level above, the Keeler TonoCare should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orientating or relocating the Keeler TonoCare .

<sup>2</sup> Over the frequency range 150kHz to 80 MHz, field strengths should be less than 10 V/m.

### 11.3 RECOMMENDED SAFE DISTANCES

#### Recommended separation distances between portable and mobile RF communications equipment and the Keeler IntelliPuff.

The Keeler TonoCare is intended for the use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Keeler TonoCare can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Keeler TonoCare as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 230MHz $d = 1.2\sqrt{p}$	80MHz to 800MHz $d = 1.2\sqrt{p}$	800MHz to 2.7GHz $d = 2.3\sqrt{p}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where  $p$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note: 1 At 80MHz and 800MHz, the higher frequency range applies.

Note 2: These guide lines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## 12. TECHNICAL SPECIFICATIONS

### TonoCare Hand Unit

<b>Hand Unit dimensions</b>	220 x 136 x 206mm (H x W x D)
<b>Hand Unit weight</b>	1.044Kg
<b>IP rating</b>	IPX0
<b>Calibrated range</b>	7mmHg to 50mmHg
<b>IOP Accuracy</b>	+/-5mmHg (95% confidence level)
<b>Working distance</b>	11mm from surface of patient's cornea to front surface of the window.
<b>Display Resolution</b>	Display resolution to 1 decimal place e.g. 12.3
<b>Display</b>	OLED 0.95"
<b>Illumination system</b>	LED, white and infrared
<b>Electric shock protection</b>	Class II (or internally powered)
<b>Complies with</b>	Electrical Safety (Medical) IEC 60601-1 IEC 60601-1-2 BS EN ISO 15004-1

The Mains plug is the means of isolating the device from the mains supply – ensure the Mains Plug is accessible at all times.

Note 1: If connecting the TonoCare to a computer, the computer shall comply with the requirements of EN 60601-1:


Note 2: Whenever the device is connected to other equipment, the combination shall comply with the requirements of EN 60601-1:

The ME Equipment includes the Hand Unit, the Docking Station, the Forehead Rest and the power supply.

**Docking Station**

<b>Docking Station dimensions</b>	153 x 155 x 183mm (H x W x D)
<b>Docking Station weight</b>	0.725Kg
<b>Electric shock protection</b>	Class II
<b>IP rating</b>	IPX0
<b>Power Supply Unit</b>	Switch mode, multi-plug type (110 – 240V) +/- 10% 350-700mA Compliant to: EN 60601-1, EN 61000-6-2, EN 61000-6-3
<b>Power supply output</b>	30 VA (12V DC 2.5A)
<b>Frequency</b>	50/60 Hz

**Environmental Conditions:**

<b>USE</b>	
	
	
Shock (without packing)	10 g, duration 6 ms
<b>STORAGE CONDITIONS</b>	
	
	
<b>TRANSPORT CONDITIONS</b>	
	
	
Vibration, sinusoidal	10 Hz to 500 Hz: 0.5g
Shock	30 g, duration 6 ms
Bump	10 g, duration 6 ms

### 13. ACCESSORIES AND SPARES

Item	Part Number
Dock	2418-P-5002
Headrest	2418-P-7000
Power supply kit	EP29-32777
TonoCare carry case	3418-P-7000
Tonometer face shield	2415-P-7038
Printer paper roll	2208-L-7008

### 14. PACKAGING AND DISPOSAL INFORMATION

#### Disposal of old electrical and electronic equipment

This symbol on the product or on its packaging and instructions indicates that this product shall not be treated as household waste.

To reduce the environmental impact of WEEE (Waste Electrical Electronic Equipment) and minimise the volume of WEEE entering landfills we encourage at product end of life that this equipment is recycled and reused.

**If you need more information on the collection reuse and recycling then please contact B2B Compliance on 01691 676124 (+44 1691 676124). (UK only).**

**Any serious incident that has occurred in relation to the device must be reported to the manufacturer and the competent authority of your Member State.**

## 15. APPENDIX

### 15.1 REPEATABILITY AND REPRODUCIBILITY

Repeatability and reproducibility of TonoCare were assessed by measuring a manometrically controlled test eye.

Repeatability was measured using a single TonoCare device and a series of approximately 50 individual readings for each of 5 pressure values spaced evenly between the 5 to 50mmHg working range. Readings were cross-referenced with a reference pressure meter and a Pulsair IntelliPuff device. Failed readings were rejected and averages of 3 subsequent readings were calculated to give a set of approximately 16 measurements at each of the 5 pressure values. Results demonstrate standard deviations ranging from 0.14mmHg to 1.11mmHg within the 5 to 50mmHg pressure range respectively.

Reproducibility was assessed by analysing measurements from three different TonoCare units by two different operators across 5 pressure values spaced evenly between the 5 to 50mmHg working range. Two measurements (an average of 4 readings) at the 5 pressure values were taken for each of the six test cases (each operator using each of the three TonoCare devices).

An Analysis of Variance (ANOVA) conducted on the data indicates a p-value of less than 0.05 and an R-square value of 98% or 99%, which signifies excellent reproducibility across operator and across devices.

### 15.2 CLINICAL PERFORMANCE DATA

#### Summary

The Keeler TonoCare Non-Contact Tonometer (NCT) was compared with the Perkins Applanation Tonometer (AT) to assess whether the TonoCare meets the requirements of ISO 8612 (comparable to ANSI Z80.10) in design compliance testing.

The Perkins AT uses the same basic principle as the Goldmann AT, namely, varying the force applied to appanate a fixed area of the cornea. Both instruments have an appanating 'cone' comprised of two prisms with apices joined together to apply an external force to the cornea to indent and flatten its surface.

There are several scientific articles referring to both instruments as reference standard tonometers and specifically the Perkins AT as the portable counterpart to the Goldmann AT (Wessels, I.F et al, 1990), (Carlos Garcia-Resua et al 2006), useful in domiciliary visits and for patients with mobility issues.

Two experienced observers acquired data from 144 qualifying eyes, measuring IOPs ranging from 7mmHg-23mmHg in 50 participants and IOPs greater than 23mmHg in 22 participants. The results of the study show that the IOP measurements taken with the TonoCare NCT when compared to the reference Perkins tonometer (AT) do not exceed the  $\pm 5$ mmHg tolerance in the three IOP ranges in 143 eyes with only 1 eye exceeding this tolerance for IOP measured  $>23$ mmHg. This falls well below the requirement that no more than 5% of the paired differences between TonoCare and the reference tonometer should be outside the  $\pm 5$ mmHg tolerance in the three IOP ranges.

Overall, the mean of IOP differences between TonoCare and Perkins AT was  $<0.01$  mmHg, with a median of  $-0.2$ mmHg, indicating that the TonoCare NCT is equivalent to the applanation tonometer.

## Methods

The study conducted was a single visit, single- centre, non-randomized, non-masked paired crossover study. The study obtained IOP measurements on each eligible eye with the TonoCare and the reference standard Perkins tonometer.

Subjects were recruited according to the following inclusion and exclusion criteria.

### Inclusion Criteria

- Subjects must be over 18 years of age
- Subjects must have healthy corneas with no contraindications for IOP measurements

### Exclusion Criteria

- Subjects with only one functional eye
- Subjects with one eye having poor or eccentric fixation
- High corneal astigmatism ( $>3D$ )
- Corneal scarring, corneal surgery (including laser corneal surgery)
- Microphthalmosis
- Buphthalmos
- Contact lens wearers
- Dry eyes
- Lid squeezers
- Nystagmus
- Keratoconus
- Any other corneal or conjunctival pathology or infection

A total of 74 eligible participants were recruited, with 2 participants (2.7%) excluded. The reason for exclusion of two participants was due to excessive blinking or anxiousness resulting in the participant holding their breath. From the included 72 participants, IOP was measured in both eyes of all participants with TonoCare and Perkins AT, giving paired IOP measurements for a total of 144 eyes.

## Results

Table 1 below gives summary of IOP characteristics of the group, showing measurements to have similar distributions.

**Table 1: Summary of TonoCare and Perkins AT IOP measurements.**

	<b>TonoCare</b>	<b>Perkins AT</b>
N, eyes (patients)	144 (72)	144 (72)
Mean IOP, mmHg	21.2	21.2
Median IOP, mmHg	18.0	17.0
SD*, mmHg	7.9	8.0
Range, mmHg	11.8 to 46.3	11.0 to 41.0
IOP 7 to 16mmHg, n (%) <sup>†</sup>	42 (29.2)	51 (35.4)
IOP 17 to 23mmHg, n (%) <sup>†</sup>	58 (40.3)	49 (34.0)
IOP >23mmHg**, n (%) <sup>†</sup>	44 (30.6)	44 (30.6)

No pairing structure is summarised in this table. \*Standard deviation.

<sup>†</sup> Only Perkins AT IOP categories are used for sub-group analyses, n is given in terms of eyes.

\*\* In order to obtain measurements in this range, an inversion procedure was performed on a subset of participants while taking IOP measurements.

Table 2 categorises the absolute differences between TonoCare and Perkins AT IOP measurements >5mmHg overall, and within 3 IOP subgroups. A difference greater than the tolerance of  $\pm 5$ mmHg occurred in 1 (0.7%) eye out of 144, well below the maximum level of 5% according to the standard.

**Table 2: Differences between TonoCare and Perkins AT IOP measurements >5mm Hg overall, and within 3 IOP subgroups.**

<b>Difference*</b>	<b>IOP Group<sup>†</sup></b>			<b>Total</b>
	<b>7 to 16mmHg</b>	<b>17 to 23mmHg</b>	<b>&gt;23mmHg**</b>	
Does not exceed $\pm 5$ mm Hg	51	49	43	143
Exceeds $\pm 5$ mmHg	0	0	1	1
Total	51	49	44	144

\* IOP TonoCare – IOP Perkins AT. <sup>†</sup>Based on Perkins AT measured IOP.

\*\* In order to obtain measurements in this range, an inversion procedure was performed on a subset of participants while taking IOP measurements.

Summary parameters of differences between pairs of TonoCare and Perkins AT IOP measurements are given in Table 3 below, for the full sample and by each IOP group. Overall the mean of IOP differences between TonoCare and Perkins AT was  $<0.01$  mmHg, with a median of  $-0.2$  mmHg. The 95% limits of agreement, based on the mean of IOP differences  $\pm 1.96 \times$  the standard deviation of the IOP differences was  $-3.4$  mmHg to  $+3.4$  mmHg.

**Table 3: Summary measures of IOP differences taken with TonoCare and Perkins AT measurements, summarised overall, and within 3 IOP subgroups.**

Summary measure [mmHg]	IOP Group <sup>†</sup>			
	7 to 16mmHg (n=51)	17 to 23mmHg (n=49)	>23mmHg** (n=44)	Overall (n=144)
Mean	0.3	0.2	-0.6	0.0
Median	0.2	0.2	-0.9	-0.2
SD*	1.3	1.5	2.1	1.7
IQR <sup>§</sup>	-0.4 to 1.2	-0.8 to 1.0	-1.8 to 0.1	-1.0 to 1.0
Range	-3 to 4	-4.0 to 4.0	-3.8 to 6.2	-4.0 to 6.2

<sup>†</sup> Based on Perkins AT measured IOP. \*Standard deviation. <sup>§</sup>Interquartile range.

\*\* In order to obtain measurements in this range, an inversion procedure was performed on a subset of participants while taking IOP measurements.

It was concluded by the investigators that there are no clinically meaningful differences in IOP measurements among the tonometers, and that the TonoCare conforms to the standard as detailed.

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Copyright © Keeler Limited 2023. Veröffentlicht in GB 2023.
Copyright © Keeler Limited 2023. Publicado en el Reino Unido 2023.
Copyright © Keeler Limited 2023. Publié au Royaume Uni en 2023.
Copyright © Keeler Limited 2023. Pubblicato nel Regno Unito 2023.
Copyright © Keeler Limited 2023. Publicado no Reino Unido 2023.
Medicel AG, Dornierstrasse 11