Pulsair desktop

Tonometer

INSTRUCTIONS FOR USE

EN


Pulsair Desktop Tonometer
Instructions For Use1

PT_BR

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CONTENTS

1. INDICATIONS FOR USE	3
1.1 BRIEF DESCRIPTION OF THE INSTRUMENT.....	3
1.2 INTENDED USE / PURPOSE OF INSTRUMENT.....	3
2. SAFETY.....	3
2.1 PHOTOTOXICITY.....	3
2.2 WARNINGS AND CAUTIONS.....	3
2.3 CONTRAINDICATION.....	5
3. CLEANING INSTRUCTIONS.....	5
3.1 CLEAN THE PUFF TUBE LENS ON A WEEKLY BASIS:	5
3.2 CLEANING TONOMETER BODY.....	5
4. POWER SUPPLY ASSEMBLY	6
4.1 SET PLUG.....	6
5. MOUNTING.....	6
6. TONOMETRY, PRESSURE VARIATIONS IN THE HUMAN EYE	6
7. NAMES OF CONTROLS AND COMPONENTS.....	7
8. SELF TEST	11
9. MEASUREMENT PROCEDURE.....	12
9.1 PREPARING THE DEVICE.....	12
9.2 PREPARING THE PATIENT.....	12
9.3 TAKING THE READING.....	13
10. PRINTING	13
10.1 USB DATA OUT.....	14
11. REPLACING THE PRINTER PAPER	14
12. USER MENU OPTIONS	14
13. CALIBRATION, MAINTENANCE AND INSPECTION.....	15
13.1 REGULAR INSPECTION.....	15
13.2 GENERAL.....	15
14. SERVICING AND CALIBRATION.....	16
15. WARRANTY	16
16. SPECIFICATIONS AND ELECTRICAL RATINGS	16
16.1 ELECTROMAGNETIC EMISSIONS.....	16
16.2 ELECTROMAGNETIC IMMUNITY.....	17
16.3 RECOMMENDED SAFE DISTANCES.....	19
17. TECHNICAL SPECIFICATIONS	19
18. ACCESSORIES AND SPARES.....	21
19. PACKAGING AND DISPOSAL INFORMATION.....	21

	Consult instructions for use		General warning sign
	Date of manufacture		Warning: Non-ionizing radiation
	Manufacturer's name and address		Warning: Electricity
	Country of manufacture		Warning: Optical radiation
	Waste Electrical and Electronic Equipment (WEEE) recycling		Warning: Floor level obstacle
	This way up		Keep dry
	Type BF applied part		Fragile
	Temperature limit		Do not use if package is damaged
UK CA 0120	United Kingdom Conformity Assessed, with the Notified Body number for SGS UK	CE 1639	Conformité Européene, with the Notified Body number for SGS Belgium NV
EC REP	Authorised representative in the European Community	CH REP	Authorised representative in Switzerland
REF	Catalogue number		Class II equipment
SN	Serial number		Atmospheric pressure limitation
MD	Medical device		Humidity limitation
	Translation		

The Keeler Pulsair Desktop Tonometer is designed and built-in conformity with Directive 93/42/EEC, Regulation (EU) 2017/745 and ISO 13485 Medical Devices Quality Management Systems.

Classification: CE / UKCA: Class IIa
FDA: Class II

The information contained within this manual must not be reproduced in whole or part without the manufacturer's prior written approval. As part of our policy for continued product development we the manufacturer reserve the right to make changes to specifications and other information contained in this document without prior notice.

This IFU is also available on the Keeler UK and Keeler USA websites.

Copyright © Keeler Limited 2023. Published in the UK 2023.

1. INDICATIONS FOR USE

These devices are intended to be used only by suitably trained and authorised healthcare professionals.

The Pulsair Desktop Tonometer should be used only by trained personnel. USA Federal law restricts this device to sale by or order of a physician.

1.1 BRIEF DESCRIPTION OF THE INSTRUMENT

It is an 'air puff' Tonometer designed to accurately measure Intra Ocular Pressure (IOP) without making contact with the surface of the eye.

Air impulse tonometry is a variant of the general applanation tonometry in which a portion of the cornea is flexed by mechanical stimuli in which the force / pressure required to produce the flexing effect is related to the intraocular pressure.

The air puff technique requires directing a calibrated quantized packet of air towards the central portion of the cornea, and the detection of the pre-defined deformation of the cornea through optical means and reflections from the corneal surface.

1.2 INTENDED USE / PURPOSE OF INSTRUMENT

The Pulsair Desktop Tonometer is indicated for measuring intraocular pressure without contacting the eye to aid in the screening and diagnosis of glaucoma.

2. SAFETY

2.1 PHOTOTOXICITY

CAUTION: The light emitted from this instrument is potentially hazardous. The longer the duration of exposure, the greater the risk of ocular damage.

While no acute optical radiation hazards have been identified for Keeler Tonometers, we recommend keeping the intensity of the light reaching the patient's retina to the minimum possible for the respective diagnosis. Children, people with aphakia and people suffering from eye conditions are most at risk. An increased risk may also occur if the retina is exposed to the same or a similar device with a visible light source within 24 hours. This applies, in particular, if the retina has been photographed with a flashbulb in advance.

Keeler Ltd shall on request, provide the user with a graph showing the relative spectral output of the instrument.

2.2 WARNINGS AND CAUTIONS

Please note that the proper and safe functioning of our instruments is only guaranteed if both the instruments and their accessories are exclusively from Keeler Ltd. The use of other

accessories may result in increased electromagnetic emissions or reduced electromagnetic immunity of the device and may lead to incorrect operation.

Observe the following precautions in order to ensure safe operation of the instruments.

WARNINGS

- Never use the instrument if visibly damaged and periodically inspect it for signs of damage or misuse.
- Check your Keeler product for signs of transport / storage damage prior to use.
- US Federal Law restricts this device to sale by or on the order of a physician or practitioner.
- The device is intended for use in various clinical settings such as hospitals, eye clinics and optometric practices.
- Only use approved Keeler power supply EP29-32777 or instrument may malfunction.
- The owner of the instrument is responsible for training personnel in its correct use.
- Never use the instrument if the ambient temperature, atmospheric pressure, and / or relative humidity are outside the limits specified in this manual.
- Do not use in the presence of flammable gases / liquids, or in an oxygen rich environment.
- This device is intended to be used only by suitably trained and authorised healthcare professionals.
- This product should not be immersed in fluid.
- The mains plug is the means of isolating the device from the mains supply. Ensure both the power switch and mains plug are always accessible.
- Do not position the equipment so that is difficult to remove the mains plug from the wall socket.

- Do not fit mains power adapter into a damaged mains outlet socket.

- Route power cords safely to eliminate risk of tripping or damage to user.

CAUTION

- Use only genuine Keeler approved parts and accessories or device safety and performance may be compromised.
- Keep out of the reach of children.
- To prevent condensation from forming, allow instrument to come to room temperature before use.
- Keeler recommends the use of disposable hygienic chinrest tissues on the chin rest before the patients place their chin on it.
- This product should be used in a room with low / dimmed lighting.

- Before using the Pulsair Desktop Tonometer, press the Clear / Demo button for 1 second to dispel any minute particles of dust or moisture which may have settled whilst the instrument was not in use.
- For indoor use only (protect from moisture).
- There are no user serviceable parts inside. Contact authorised service representative for further information.
- Follow guidance on cleaning / routine maintenance to prevent personal injury / damage to equipment.
- Failure to carry out recommended routine maintenance as per the instructions in this IFU may reduce the operational lifetime of the product.
- At product end of life dispose of in accordance with local environmental guidelines (WEEE).

2.3 CONTRAINDICATION

There is no restriction to patient population this device can be used with other than those outlined in the contraindications stated below.

Accuracy of IOP measurements is known to be affected by variations and changes in corneal rigidity due to differences in corneal thickness, intrinsic structural factors or corneal refractive surgery. It is recommended that these factors are considered during IOP measurement.

3. CLEANING INSTRUCTIONS

3.1 CLEAN THE PUFF TUBE LENS ON A WEEKLY BASIS:

1. Moisten a cotton bud with Isopropyl Alcohol.
2. Move the tip of the bud around the lens in a circular motion.
3. After one circle the bud should be discarded to avoid smearing on the lens.
4. Look at the Puff Tube lens from the patient's side, if traces of tear film can still be seen, repeat above steps until clear.

Note: Care should be taken not to damage the Puff Tube assembly during cleaning.

CAUTION: Never use a dry cotton bud or tissue to clean the Puff Tube lens. Never use a silicone impregnated cloth or tissue to clean the Puff Tube lens.

3.2 CLEANING TONOMETER BODY

Only manual non-immersion cleaning as described should be used for this tonometer. Do not use corrosive products. Do not autoclave or immerse in cleaning fluids. Always disconnect power supply from source before cleaning.

1. Wipe the external surface with a clean absorbent, non-shedding cloth dampened with de-ionised water / detergent solution (2% detergent by volume) or water / isopropyl alcohol solution (70% IPA by volume). Avoid optical surfaces.
2. Ensure that excess solution does not enter the instrument. Use caution to ensure cloth is not saturated with solution.
3. Surfaces must be carefully hand-dried using a clean non-shedding cloth.
4. Safely dispose of used cleaning materials.

4. POWER SUPPLY ASSEMBLY

4.1 SET PLUG

Replace the blanking plate with the appropriate mains plug adapter if required, or use IEC 60320 TYPE 7 connector (not supplied).

5. MOUNTING

Your Pulsair Desktop Tonometer is designed to be used on a sturdy flat surface, for example, a medical instrument table or a purpose designed refraction system tabletop.

Choose carefully the intended location for your Pulsair Desktop Tonometer with particular consideration to health and safety aspects, for example, the routing of the power lead, and its position in regard to the user and the patient.

Your Pulsair Desktop Tonometer has four anti slip rubber feet, check these are located well within the edge of your intended flat surface to ensure there is no possibility of the Pulsair Desktop Tonometer becoming dislodged and causing injury to a user or patient.

Your Pulsair Desktop Tonometer has an adjustable chinrest, however, we recommend that to ensure maximum patient comfort you use it on a height adjustable table that allows wheelchair access for any wheelchair bound patients.

6. TONOMETRY, PRESSURE VARIATIONS IN THE HUMAN EYE

The Keeler Pulsair Desktop Tonometer measures intra-ocular pressure by automatically releasing a gentle puff of air onto the cornea. This is known as an event.

A single reading can sometimes be misleading as the IOP will vary because of pulse, respiratory and diurnal fluctuations. In addition, blinking, squeezing, fluid intake, physical activity, body position and even the direction of gaze can influence IOP.

Up to 4 readings may be required to reduce the impact of these variants to obtain a constant IOP.

Pulsair Desktop Tonometer software will recognise the readings and emit a sound notification when two consecutive readings are ± 1 mmHg of each other indicating that further measurements may not be required.

7. NAMES OF CONTROLS AND COMPONENTS

1 Travelling Lock

The Pulsair Desktop Tonometer is fitted with a Travelling Lock to protect the moving parts from damage when in transit or when the tonometer may be subject to sudden movement or shock.


To release the Travelling Lock rotate it counter-clockwise until the Joystick Assembly is free to move.

To lock the moving section in place, align the Moving Platform with the Base Unit and carefully screw the Travelling Lock clockwise until it is firmly located.

Note: do not over tighten the Travelling Lock.


2 Chinrest and Chinrest Height Adjuster

Rotate the Chinrest Adjuster to increase or decrease the height of the Chinrest until the patient's outer canthus is in line with the indicator line on the Chinrest vertical bar.

3 On / Off Push Button

Press the On / Off button to switch on the tonometer and ready for use. The display will progress quickly through the following stages during system initialisation.

When the system is ready to use, the display will read 'L: READY' or 'R: READY' depending on whether the tonometer is in position to test the left eye or the right eye.

'R' indicates it is in the right eye position – the display could also read 'L' for the left eye position.

4 Start / Stop Button

Pressing the Start / Stop button while the Pulsair Desktop Tonometer is running will stop the pump and put the system into Standby Mode, the display will show 'STANDBY'. The Start / Stop button will reset the memory buffer, clearing all readings stored. Pressing the Start / Stop button while the unit is in Standby Mode, the Pulsair Desktop Tonometer will start with the pump running and the system is initialised ready for use. The tonometer will go through all the displays.

5 Motion Sensitive Wakeup Mechanism

The Pulsair Desktop Tonometer is equipped with a motion sensor. Upon moving the Moving Platform from left to right or vice versa, the motion sensitive switch is activated to start the pump and the rest of the system from Standby Mode so that it is ready for measurement use. The tonometer will go through all the displays. After two minutes of inactivity the system automatically switches off to save energy.

6 Clear / Demo Button

Pressing the Clear / Demo button momentarily clears all records of previous readings, and the tonometer is set to its default setting. If Clear / Demo button is pressed and held for more than a second, the tonometer initialises and fires a demo puff to demonstrate to the patient the softness of the air puff.

7 Print / Menu Button

Pressing the Print / Menu button momentarily will print the results on the integrated printer and send the data to the USB data port. If Print / Menu button is pressed and held for more than a second, the User


Menu Options is entered.

8 Easy Pulse Button

Pressing and holding the Easy Pulse button for over a second activates the Easy Pulse Mode.

This enables the tonometer to fire when it is difficult to get a measurement, for example, with a damaged or scarred cornea. This action is accompanied by sound and an additional momentary indication on the display. On pressing any other button or changing the eye, the tonometer resets to normal mode.

9 USB Transfer Port

This is used for transferring IOP data to an external device, such as a personal computer, and for use by a qualified technician when calibrating the Pulsair Desktop Tonometer or performing a software revision.

10 Power Input Connection / Disconnection

To connect, insert the low voltage power lead into the power input socket. To disconnect, remove the low voltage power lead from the power input socket.

11 Joystick

The rotating Joystick controls the elevation of the Measurement Head and back-front movements of the Moving Platform of the Pulsair Desktop Tonometer.

12 Printer Cover

Access to the printer paper is via the Printer Cover, pull the lip on the top of the cover and gently pull towards you to open the Printer Cover.

13 Alignment Screen

The Alignment Screen allows the user to visualise the patient's eyes to correctly align the tonometer with the centre of the cornea prior to measurement. Alignment is performed using the Moving Platform of the Pulsair Desktop Tonometer and the Joystick for final alignment.

The Pulsair Desktop Tonometer will fire automatically when correctly positioned and aligned. See section 8 for detailed instructions on the alignment process.

14 Alignment Screen Adjustment Control Wheels

Colour adjustment

Brightness adjustment


Contrast adjustment

15 Measurement Display Screen

This screen will display the eye measured, the average reading and individual readings from either the left (L) and right (R) eyes. Left and right detection is automatic.

After the first reading is taken the display shows the measured

IOP. After each of the consecutive readings are taken the display shows the individual reading and the average of up to the last four readings per eye.

In addition to IOP readings, the Pulsair Desktop Tonometer also displays several messages on the display when a measurement is not detected for a number of reasons. In such cases, the display may read as follows:

- a) < 5: This is shown when the reading is lower than 5mmHg, in which case no valid reading is displayed. Signified by long sound.

- b) >25: This is shown when applanation is detected with soft-puff and the reading is greater than 25mmHg, in this case the tonometer displays '>25' and switches to hard puff. The tonometer shall stay in hard puff mode until two consecutive readings are less than 20mmHg.

- c) >50: This is shown when the reading is greater than 50mmHg, in which case no valid reading is displayed. Signified by long sound.

Note: The displayed average is rounded to the nearest whole number or displayed to one decimal place depending on the user setting accessed via the User Menu Options.

The displayed running average is based on the actual readings which are taken to one decimal place. For example, readings of 15.4, 16.3, 14.2 and 16.9 are averaged by adding them together, which equals 62.8, and dividing by the number of readings taken, 4. This gives a final figure of 15.7, or 16 depending on user settings. Please note independent readings are displayed as whole numbers.

When all the required readings have been taken the figure displayed is the IOP that is recorded for the patient. When two consecutive readings are within 1mmHg an audible sound will be emitted indicating that sufficient readings may have been taken.

8. SELF TEST

A self test program can be run for diagnostic purpose by selecting the Menu Option. The display will show 'RUN SELF TEST'. By pressing the Clear / Demo button self test is entered and the display will show 'SELF TEST RUNNING', and this will run for 45 seconds. At the end of the test the resulting data is sent to the printer and the display is cleared.

9. MEASUREMENT PROCEDURE

9.1 PREPARING THE DEVICE

1. Plug in the Power Supply Cord to the tonometer. The Power Socket is located on the right-hand side of the tonometer.
2. Turn the tonometer on using the On / Off push switch located on the front of the tonometer. The tonometer will initialise ready for use.
3. Remove the protective dust cap from the Puff Tube.
4. Unscrew the Travelling Lock if secured.
5. Using the Joystick bring the moving part of the Pulsair Desktop Tonometer back towards you and across to the left (to measure the right eye first).
6. Before using the Pulsair Desktop Tonometer, press the Clear / Demo button for 1 second to dispel any minute particles of dust or moisture which may have settled on the tonometer whilst not in use.


9.2 PREPARING THE PATIENT

Before using the Pulsair Desktop Tonometer you should make your patient feels at ease and ensure they are in an optimum reading location, preferably with their head supported. This is because apprehension and nervousness may adversely affect the readings obtained. Follow the points outlined below to achieve this:

1. Ask the patient to remove their contact lenses or spectacles if worn and to blink and breathe normally.
2. Ensure that the patient is comfortable and in a relaxed position.
3. To reassure the patient, you can demonstrate the procedure, using the Clear / Demo button, on the back of the patient's hand prior to taking a reading.
4. Place disposable chinrest tissue on the Chinrest.
Ask the patient to place their head on the Chinrest.

5. Adjust the Chinrest height so that the outer canthus is aligned with the marker on the Chinrest vertical support.


Before taking a reading you should:

1. Ask the patient to blink to ensure a good and reflective tear film.
2. Ensure the patient and tonometer optics are not positioned under direct lighting (i.e. spot lights or sunlight).
3. Ensure the patient's eyes are fully opened. This helps to prevent squeezing, where the patient unconsciously tenses their eyelids and increases IOP.
4. Throughout the reading process, you should allow the patient to blink at intervals to maintain the corneal tear film.

9.3 TAKING THE READING

1. Holding the Joystick, move the Moving Platform with the other hand until the patients eye to be measured is visible and central on the Alignment Screen. Height adjustment is obtained by rotating the Joystick. If there is insufficient travel re-check the patients head height in the Chinrest and try again using the Joystick
2. Carefully move the tonometer towards the patient until the external eye image becomes the cross hairs alignment target.
3. Using the Joystick, focus the cross-hair alignment target until the correct focused position triggers the tonometer to fire automatically.

Note: a red hue may be visible on the screen, this does not affect the function of the device.

4. Remain in the firing position until the Pulsair Desktop Tonometer stops reading after four readings per eye are taken. The Pulsair Desktop Tonometer will emit a short beep when sufficient readings have been taken i.e. when two successive readings are within +/- 1mmHg of each other.
5. If a reading is recorded as a non-event or bad event, a long tone will be heard.
6. The first reading will be the measured value; successive readings will display the running average IOP. Outlying or spurious readings will automatically be excluded from the calculation.

10. PRINTING

The results can be printed by pressing the print button at any time. Printing the results does not clear the print memory buffer.


10.1 USB DATA OUT

Pressing the print button sends a flat file to the USB port in the following format:

[R] = xx, xx, xx, xx, [xx.x]


[L] = xx, xx, xx, xx, [xx.x]

This data may be able to be imported into other applications, for details of how to perform this please consult your third party program support team.

11. REPLACING THE PRINTER PAPER

1. Access to the printer paper is via the Printer Cover, pull the lip on the top of the Cover and gently pull towards you to open the Printer Cover. Remove the empty paper roll.
2. Place the new roll of paper into the paper holder, making sure the free end is loose at the top of the roll, otherwise it won't print.
3. Feed the free end of paper through the gap in the Cover.
4. Close the Cover.

12. USER MENU OPTIONS

1. With the tonometer switched on, press and hold the Print / Menu button for more than 3 seconds to enter the User Menu.
2. The display will show the first Menu Option (Buzzer Control) and the current selection (i.e.. [BUZZER ON] or [BUZZER OFF])

- Pressing the Clear / Demo button for more than 1 second, cycles through the user Change Options (shown in table.)
- Pressing the Print / Menu button for more than 1 second moves you forward to the next Menu Option (shown in table), in this case the Desktop Level.
- Use the Clear / Demo button to make your preferred selection.
- Continue to repeat steps 4 and 5 until 'OK' is displayed. Your Pulsair Desktop Tonometer is now ready to use with your preferred settings.

Menu Option	Display	Change Options
Buzzer Control	BUZZER ON	ON / OFF
IOP Format	DISPLAY XX	XX / XX.X
Display Screen	CONTRAST 0	0 - 20
Full Self Test	RUN SELF TEST?	RUNNING SELF TEST

13. CALIBRATION, MAINTENANCE AND INSPECTION

Keeler recommends this routine maintenance be carried out by the user frequently to ensure safe and accurate measurement. In the event of the device being outside of the calibration tolerances, it is important to send the device back to Keeler Ltd. or your local dealer for repair and re-calibration.

13.1 REGULAR INSPECTION

Inspect your power supply unit and cable for damage regularly.

Before inspecting, disconnect the power supply from the Pulsair Desktop Tonometer and the mains.

If the outer insulation of the cable appears to be damaged discontinue use immediately. Contact your local dealer for a replacement.

13.2 GENERAL

Keep the tonometer free from dust.

If the Pulsair Desktop Tonometer is to remain unused for any length of time, press the On / Off Push button switch to 'Off' and remove the power supply. Use the dust cover to protect the tonometer.

14. SERVICING AND CALIBRATION

Keeler recommends an annual calibration for the Tonometer. Do not modify this equipment without authorisation of the manufacturer.

This must be performed by an authorised Pulsair service centre or distributor. The unit performs a self function check when switched on and will indicate if a fault is found.

There are no user serviceable parts in this instrument. Service manuals will be available to authorised Keeler service centres and Keeler trained service personnel.

15. WARRANTY

Your Keeler product is guaranteed for 2 years and will be replaced, or repaired free of charge subject to the following:

- Any fault due to faulty manufacture.
- The instrument and accessories have been used in compliance with these instructions.
- Proof of purchase accompanies any claim.

The manufacturer declines any and all responsibility and warranty coverage should the instrument be tampered with in any manner or should routine maintenance be omitted or performed in manners not in accordance with these manufacturer's instructions.

There are no user serviceable parts in this instrument. Any servicing or repairs should only be carried out by Keeler Ltd. or by suitably trained and authorised distributors. Service manuals will be available to authorised Keeler service centres and Keeler trained service personnel.

16. SPECIFICATIONS AND ELECTRICAL RATINGS

The Keeler Pulsair Desktop Tonometer is a medical electrical instrument. The instrument requires special care concerning electromagnetic compatibility (EMC). This Section describes its suitability in terms of electromagnetic compatibility of this instrument. When installing or using this instrument, please read carefully and observe what is described here.

Portable or mobile-type radio frequency communication units may have an adverse effect on this instruments, resulting in malfunctioning.

16.1 ELECTROMAGNETIC EMISSIONS

Guidance and manufacturer's declaration – electromagnetic emissions

The Keeler Pulsair Desktop tonometer is intended for use in the electromagnetic environment specified below. The customer or user should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Keeler Pulsair Desktop Tonometer uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Keeler Pulsair Desktop Tonometer is suitable for use in all establishments, other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2		
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

16.2 ELECTROMAGNETIC IMMUNITY

Guidance and manufacturer's declaration – electromagnetic immunity

The Keeler Pulsair Desktop Tonometer is intended for use in the electromagnetic environment specified below. The customer or user should assure that it is used in such an environment.

Immunity test	IEC 55015 Test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD). IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst. IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge. IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines. IEC 61000-4-11	$U_T = 0\%$ 0.5 cycle (0, 45, 90, 135, 180, 225, 270, 315°) $U_T = 0\%$; 1 cycle $U_T = 70\%$; 25/30 cycles (@ 0°) $U_T = 0\%$; 250/300 cycle	$U_T = 0\%$ 0.5 cycle (0, 45, 90, 135, 180, 225, 270, 315°) $U_T = 0\%$; 1 cycle $U_T = 70\%$; 25/30 cycles (@ 0°) $U_T = 0\%$; 250/300 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Keeler Pulsair Desktop Tonometer requires continued operation during power mains interruptions, it is recommended that the charger be powered from an uninterruptible power supply.

Immunity test	IEC 55015 Test level	Compliance level	Electromagnetic environment – guidance
Power frequency (50/60 Hz) Magnetic field. IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at a level characteristic of a typical location in a typical professional healthcare facility environment.

Note: U_T is the a. c. mains voltage prior to application of the test level.

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment – guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the Keeler Pulsair Desktop Tonometer, including cables, than the recommended separation distances calculated from the equation applicable to the frequency of the transmitter.
		Recommended separation distance	
Conducted RF IEC 61000-4-6	6 Vrms 1	6 V	$d = 1.2 \sqrt{p}$
Radiated RF IEC 61000-4-3	10 V/m 80MHz to 2.7GHz	10 V/m	$d = 1.2 \sqrt{p}$ 80MHz to 800 MHz $d = 2.3 \sqrt{p}$ 800MHz to 2.7GHz
			Where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ¹ , should be less than the compliance level in each frequency range. ²
			 Interference may occur in the vicinity of equipment marked with this symbol.

Note 1: At 80MHz and 800MHz, the higher frequency range applies.

Note 2: These guide lines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people

¹ Field strengths from fixed transmitters, such as base stations (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Keeler Pulsair Desktop Tonometer is used exceeds the applicable RF compliance level above, the Keeler Pulsair Desktop Tonometer should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as

re-orientating or relocating the Keeler Pulsair Desktop Tonometer.

² Over the frequency range 150kHz to 80 MHz, field strengths should be less than 10 V/m.

16.3 RECOMMENDED SAFE DISTANCES

Recommended separation distances between portable and mobile RF communications equipment and the Keeler Desktop.

The Keeler Pulsair Desktop Tonometer is intended for the use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Keeler Pulsair Desktop Tonometer can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Keeler Pulsair Desktop Tonometer as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 230MHz $d = 1.2\sqrt{p}$	80MHz to 800MHz $d = 1.2\sqrt{p}$	800MHz to 2.7GHz $d = 2.3\sqrt{p}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note: 1 At 80MHz and 800MHz, the higher frequency range applies.

Note 2: These guide lines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.





17. TECHNICAL SPECIFICATIONS



Dimensions	450 x 435 x 245mm (H x D x W)
Weight	16kg
Calibrated range	5mmHg to 50mmHg
Repeatability (Average coefficient of variation)	<5%
Accuracy	+/-5mmHg (95% confidence level)*

Working distance	20mm from surface of patient's cornea to front surface of first lens. This equates to a nominal distance of 15mm from the front of the puff tube shroud to the front surface of the patient's cornea
Displayed scale	Single line 16 character alphanumeric display
Illumination system	LED infra-red
Complies with	Electrical Safety (Medical) IEC 60601-1 Electromagnetic Compatibility IEC 60601-1-2, BS EN ISO 15004-1, BS EN ISO 15004-2
Power Supply Unit	Switch mode, (110-240V)+/- 10% multi plug type compliant to EN 60601-1, EN 61000-6-2, EN 61000-6-3
Power supply output	30 VA (12V DC 2.5A)
Frequency	50/60 Hz

*In house trial carried out on 20 participants

Environmental Conditions:

USE		
		
Shock (without packing)	10 g, duration 6 ms	
STORAGE CONDITIONS		
		
TRANSPORT CONDITIONS		

	  
Vibration, sinusoidal	10 Hz to 500 Hz: 0.5g
Shock	30 g, duration 6 ms
Bump	10 g, duration 6 ms

18. ACCESSORIES AND SPARES

Item	Part Number
Chinrest paper	3104-L-8201
Printer paper	2208-L-7008
Pulsair Desktop Tonometer dust cover	EP39-70435

19. PACKAGING AND DISPOSAL INFORMATION

Disposal of old electrical and electronic equipment

This symbol on the product or on its packaging and instructions indicates that this product shall not be treated as household waste.

To reduce the environmental impact of WEEE (Waste Electrical Electronic Equipment) and minimise the volume of WEEE entering landfills we encourage at product end of life that this equipment is recycled and reused.

If you need more information on the collection reuse and recycling then please contact B2B Compliance on 01691 676124 (+44 1691 676124). (UK only).

Any serious incident that has occurred in relation to the device must be reported to the manufacturer and the competent authority of your Member State.

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Copyright © Keeler Limited 2023. Publicado no Reino Unido em 2023.
Medicel AG, Dornierstrasse 11
SGS United Kingdom Ltd