**Pulsair** **intelliPuff**

## Non-Contact Tonometer

INSTRUCTIONS FOR USE

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## CONTENTS

<b>1. INDICATIONS FOR USE</b> .....	<b>3</b>
1.1 BRIEF DESCRIPTION OF THE INSTRUMENT .....	3
1.2 INTENDED USE / PURPOSE OF INSTRUMENT .....	3
<b>2. SAFETY</b> .....	<b>3</b>
2.1 PHOTOTOXICITY .....	3
2.2 WARNINGS AND CAUTIONS .....	4
2.3 CONTRAINDICATION .....	5
<b>3. CLEANING INSTRUCTIONS</b> .....	<b>5</b>
3.1 CLEAN THE PUFF TUBE LENS ON A WEEKLY BASIS: .....	5
3.2 CLEANING TONOMETER BODY .....	6
<b>4. POWER SUPPLY ASSEMBLY</b> .....	<b>6</b>
4.1 SET PLUG .....	6
<b>5. WALL MOUNTING</b> .....	<b>6</b>
<b>6. TONOMETRY, PRESSURE VARIATIONS IN THE HUMAN EYE</b> .....	<b>7</b>
<b>7. NAMES OF CONTROLS AND COMPONENTS</b> .....	<b>7</b>
<b>8. MEASUREMENT PROCEDURE</b> .....	<b>10</b>
8.1 PREPARING THE DEVICE .....	10
8.2 PREPARING THE PATIENT .....	11
8.3 TAKING THE READING .....	11
<b>9. DISPLAY EXAMPLES</b> .....	<b>13</b>
<b>10. PRINTING</b> .....	<b>14</b>
10.1 SAMPLE PRINT .....	14
<b>11. REPLACING THE PRINTER PAPER</b> .....	<b>14</b>
<b>12. USER MENU OPTIONS</b> .....	<b>14</b>
<b>13. CALIBRATION, MAINTENANCE AND INSPECTION</b> .....	<b>15</b>
13.1 REGULAR INSPECTION .....	15
13.2 GENERAL .....	16
<b>14. SERVICING AND CALIBRATION</b> .....	<b>16</b>
<b>15. WARRANTY</b> .....	<b>16</b>
<b>16. SPECIFICATIONS AND ELECTRICAL RATINGS</b> .....	<b>16</b>
16.1 ELECTROMAGNETIC EMISSIONS .....	17
16.2 ELECTROMAGNETIC IMMUNITY .....	17
16.3 RECOMMENDED SAFE DISTANCES .....	19
<b>17. TECHNICAL SPECIFICATIONS</b> .....	<b>20</b>
<b>18. ACCESSORIES AND SPARES</b> .....	<b>21</b>
<b>19. PACKAGING AND DISPOSAL INFORMATION</b> .....	<b>21</b>

	Consult instructions for use		General warning sign
	Date of manufacture		Warning: Non-ionizing radiation
	Manufacturer's name and address		Warning: Electricity
	Country of manufacture		Warning: Optical radiation
	Waste Electrical and Electronic Equipment (WEEE) recycling		Warning: Floor level obstacle
	This way up		Keep dry
	Type BF applied part		Fragile
	Temperature limit		Do not use if package is damaged
<b>UK CA 0120</b>	United Kingdom Conformity Assessed, with the Notified Body number for SGS UK	<b>CE1639</b>	Conformité Européenne, with the Notified Body number for SGS Belgium NV
<b>EC</b>   <b>REP</b>	Authorised representative in the European Community	<b>CH</b>   <b>REP</b>	Authorised representative in Switzerland
<b>REF</b>	Catalogue number		Class II equipment
<b>SN</b>	Serial number		Atmospheric pressure limitation
<b>MD</b>	Medical device		Humidity limitation
	Translation		

The Keeler Pulsair Intellipuff Tonometer is designed and built-in conformity with Directive 93/42/EEC, Regulation (EU) 2017/745 and ISO 13485 Medical Devices Quality Management Systems.

Classification: CE / UKCA: Class IIa  
FDA: Class II

The information contained within this manual must not be reproduced in whole or part without the manufacturer's prior written approval. As part of our policy for continued product development we the manufacturer reserve the right to make changes to specifications and other information contained in this document without prior notice.

This IFU is also available on the Keeler UK and Keeler USA websites.

Copyright © Keeler Limited 2023. Published in the UK 2023.

## 1. INDICATIONS FOR USE

These devices are intended to be used only by suitably trained and authorised healthcare professionals.

**The Pulsair IntelliPuff Non-Contact Tonometer should be used only by trained personnel. USA Federal law restricts this device to sale by or order of a physician.**

### 1.1 BRIEF DESCRIPTION OF THE INSTRUMENT

It is an 'air puff' Tonometer designed to accurately measure Intra Ocular Pressure (IOP) without making contact with the surface of the eye.

This device is intended to be used by a trained healthcare professional only in a healthcare environment. Air impulse tonometry is a variant of the general applanation tonometry in which a portion of the cornea is flexed by mechanical stimuli in which the force / pressure required to produce the flexing effect is related to the intraocular pressure.

The air puff technique requires directing a calibrated quantized packet of air towards the central portion of the cornea, and the detection of the pre-defined deformation of the cornea through optical means and reflections from the corneal surface.

### 1.2 INTENDED USE / PURPOSE OF INSTRUMENT

The Pulsair IntelliPuff Non-Contact Tonometer is indicated for measuring intraocular pressure without contacting the eye to aid in the screening and diagnosis of glaucoma.

## 2. SAFETY

### 2.1 PHOTOTOXICITY

**CAUTION: The light emitted from this instrument is potentially hazardous. The longer the duration of exposure, the greater the risk of ocular damage.**

**While no acute optical radiation hazards have been identified for Keeler Tonometers, we recommend keeping the intensity of the light reaching the patient's retina to the minimum possible for the respective diagnosis. Children, people with aphakia and people suffering from eye conditions are most at risk. An increased risk may also occur if the retina is exposed to the same or a similar device with a visible light source within 24 hours. This applies, in particular, if the retina has been photographed with a flashbulb in advance.**

**Keeler Ltd shall on request, provide the user with a graph showing the relative spectral output of the instrument.**

### 2.2 WARNINGS AND CAUTIONS

Please note that the proper and safe functioning of our instruments is only guaranteed if both the instruments and their accessories are exclusively from Keeler Ltd. The use of other accessories may result in increased electromagnetic emissions or reduced electromagnetic immunity of the device and may lead to incorrect operation.

Observe the following precautions in order to ensure safe operation of the instruments.

#### WARNINGS

- Never use the instrument if visibly damaged and periodically inspect it for signs of damage or misuse.
- Check your Keeler product for signs of transport / storage damage prior to use.
- US Federal Law restricts this device to sale by or on the order of a physician or practitioner.
- The device is intended for use in various clinical settings such as hospitals, eye clinics and optometric practices.
- Only use approved Keeler power supply EP29-32777 or instrument may malfunction.
- Do not remove the labels covering the key holes unless wall mounting the Pulsair IntelliPuff.
- The owner of the instrument is responsible for training personnel in its correct use.
- Never use the instrument if the ambient temperature, atmospheric pressure, and / or relative humidity are outside the limits specified in this manual.
- Do not use in the presence of flammable gases / liquids, or in an oxygen rich environment.
- This device is intended to be used only by suitably trained and authorised healthcare professionals.
- This product should not be immersed in fluid.
- The mains plug is the means of isolating the device from the mains supply. Ensure both the power switch and mains plug are always accessible.
- Do not position the equipment so that is difficult to remove the mains plug from the wall socket.

- Do not fit mains power adapter into a damaged mains outlet socket.

- Route power cords safely to eliminate risk of tripping or damage to user.

#### CAUTION

- Use only genuine Keeler approved parts and accessories or device safety and performance may be compromised.
- Keep out of the reach of children.

- To prevent condensation from forming, allow instrument to come to room temperature before use.
- Only mount on wall according to Keeler Instructions.
- This product should be used in a room with low / dimmed lighting.
- Before using the Pulsair IntelliPuff Tonometer, press the Demo button for 1 second to dispel any minute particles of dust or moisture which may have settled whilst the instrument was not in use.
- For indoor use only (protect from moisture).
- There are no user serviceable parts inside. Contact authorised service representative for further information.
- Follow guidance on cleaning / routine maintenance to prevent personal injury / damage to equipment.
- Failure to carry out recommended routine maintenance as per the instructions in this IFU may reduce the operational lifetime of the product.
- At product end of life dispose of in accordance with local environmental guidelines (WEEE).

### 2.3 CONTRAINDICATION

There is no restriction to patient population this device can be used with other than those outlined in the contraindications stated below.

Accuracy of IOP measurements is known to be affected by variations and changes in corneal rigidity due to differences in corneal thickness, intrinsic structural factors or corneal refractive surgery. It is recommended that these factors are considered during IOP measurement.

## 3. CLEANING INSTRUCTIONS

### 3.1 CLEAN THE PUFF TUBE LENS ON A WEEKLY BASIS:

1. Moisten a cotton bud with Isopropyl Alcohol.
2. Move the tip of the bud around the lens in a circular motion.
3. After one circle the bud should be discarded to avoid smearing on the lens.
4. Look at the Puff Tube lens from the patient's side, if traces of tear film can still be seen, repeat above steps until clear.

**Note:** Care should be taken not to damage the Puff Tube assembly during cleaning.

**CAUTION:** Never use a dry cotton bud or tissue to clean the Puff Tube lens. Never use a silicone impregnated cloth or tissue to clean the Puff Tube lens.

### 3.2 CLEANING TONOMETER BODY

Only manual non-immersion cleaning as described should be used for this tonometer. Do not autoclave or immerse in cleaning fluids. Always disconnect power supply from source before cleaning.

1. Wipe the external surface with a clean absorbent, non-shedding cloth dampened with de-ionised water / detergent solution (2% detergent by volume) or water / isopropyl alcohol solution (70% IPA by volume). Avoid optical surfaces.
2. Ensure that excess solution does not enter the instrument. Use caution to ensure cloth is not saturated with solution.
3. Surfaces must be carefully hand-dried using a clean non-shedding cloth.
4. Safely dispose of used cleaning materials.

## 4. POWER SUPPLY ASSEMBLY

### 4.1 SET PLUG

Replace the blanking plate with the appropriate mains plug adapter if required, or use IEC 60320 TYPE 7 connector (not supplied).

## 5. WALL MOUNTING


Your Pulsair IntelliPuff is supplied with a sturdy wall mounting bracket.

1. The bracket has four holes allowing it to be securely fixed to an appropriate wall or vertical surface.
2. Choose carefully the intended location for your IntelliPuff with particular consideration to health and safety aspects, for example, the routing of the power lead, and its position in regard to the user and the patient.
3. Use the wall mounting bracket as a template and clearly mark the position of the holes in the wall. Ensure that there are no live utilities where you are to drill.
4. Drill the appropriate size holes for the screws and rawplugs supplied.
5. Attach the plate securely to the wall.
6. Remove labels covering the key holes.
7. Carefully locate the mounting pegs of the plate to the key holes on the back of your Pulsair IntelliPuff and lower into its final and secure position.

## 6. TONOMETRY, PRESSURE VARIATIONS IN THE HUMAN EYE

The Keeler Pulsair IntelliPuff Tonometer measures intra-ocular pressure by automatically releasing a gentle puff of air onto the cornea. This is known as an event.

A single reading can sometimes be misleading as the IOP will vary because of pulse, respiratory and diurnal fluctuations. In addition, blinking, squeezing, fluid intake, physical activity, body position and even the direction of gaze can influence IOP.

Up to 4 readings may be required to reduce the impact of these variants to obtain a constant IOP.

Pulsair IntelliPuff Tonometer software will recognise the readings and emit a sound notification when two consecutive readings are  $\pm 1$  mmHg of each other indicating that further measurements may not be required.

## 7. NAMES OF CONTROLS AND COMPONENTS

### 1 On / Off Push Button

To turn the Pulsair IntelliPuff on or off, push the On / Off button – a green LED will indicate the unit is on.

### 2 Test Eye

This is useful for user training; it will not return an IOP measurement.

### 3 Printer Active LED

When lit, this indicates the printer is activated; replacing the handset in the cradle will start the print. Alternatively, a print can be made at any time using the print button on the handset.

### 4 Printer Cover

Access to the printer paper is via the Printer Cover, pull the lip on the top of the cover and gently pull towards you to open the Printer Cover.

### 5 Serial Port

The Serial Port is used for calibration, systems checking and data output (located on rear of instrument).

### 6 Power Input


Insert the low voltage power lead (located on rear of instrument) using only Keeler power supplies.

### 7 Forehead Rest

Push to release or push to return the Forehead Rest to its discrete position.

### 8 Right (OD) / Left (OS) Indicators

These will indicate the eye to be measured; the OD / OS button will toggle between these indicators.

## 9 Display

The display shows the recorded IOP reading and the averaged IOP reading.

After the first reading is taken the display shows the measured IOP. After each of the consecutive reading is taken the display shows the average of the readings taken so far, i.e. the first figure displayed is the actual reading, the second figure is an average of the first two readings etc., up to a maximum of 4 readings per eye.

**Note:** The displayed figure is rounded to the nearest whole number or displayed to one decimal place depending on the user setting accessed via the User Menu Options.

The displayed average is based on the readings which are taken to one decimal place. For example, readings of 15.4, 16.3, 14.2 and 16.9 are averaged by adding them together which equals 62.8 and dividing by the number of readings taken, 4. This gives a final figure of 15.7 or 16 depending on user settings.

When all the required readings have been taken the figure displayed is the IOP that is recorded for the patient. When two consecutive readings are within 1mmHg an audible sound will be heard indicating that sufficient readings have been taken.

## 10 Print / Menu button

A press of less than one second will print the acquired data; press and hold for more than 3 seconds to access the User Menu Options. Refer to page 14 for full instructions on the User Menu Options.

## 11 Review Button / Easy Pulse Button

The 'R' Review button is dual function:

- Review – It allows the reader to review readings taken.
- Easy Pulse Mode – In the event of difficulties in firing. For example with a damaged or scarred cornea, it will override the firing parameters to ease taking the measurement.

Press the Review button. The display shows the readings taken in the order they were taken, the final figure displayed is the cumulative average, the IOP.

The Pulsair memory retains a rolling four readings, per eye. New readings automatically replace the oldest.

To review the other eye, press the OD / OS button once and then press the Review button.

To clear the memory you can either replace the Hand Unit in the holster and remove again or press the Demo button.

To initiate Easy Pulse Mode hold the Review button for greater than one second; the display will show 'easy', it will beep once and the Pulsair Intellipuff tonometer will be ready to use on the difficult eye. Pressing any button, returning the Hand Unit to the Cradle, or performing a manual reset by pressing the button in the Cradle Well, will return the Pulsair IntelliPuff to its previous settings.

## 12 Demo Button

To reassure the patient, you can demonstrate the procedure, using the Demo button, on the back of the patient's hand prior to taking a reading.

## 13 OD / OS Button 'Menu Change Button'

This toggles between recording data for the left or right eyes. This button is also used to toggle through the User Menu Options when in Menu Mode, refer to page 14 for full instructions on the User Menu Options.

## 14 Eyepiece

The Eyepiece allows the user to view the patient's eye and align the targeting system.

## 15 Puff Tube and Lens

The Puff Tube and Puff Lens are the parts of the Pulsair Intellipuff through which the Pulsair Intellipuff is aligned and a gentle puff of air is emitted.

## 16 Alignment LEDs


The two green LEDs located on the front of the Hand Unit act as a guide when you are lining up the patient's eye to take a reading.

## 8. MEASUREMENT PROCEDURE

### 8.1 PREPARING THE DEVICE

1. Plug in the Power Supply Cord to the tonometer. The Power Socket is located at the rear of the tonometer.
2. Turn the tonometer on using the On / Off push switch located on the front of the tonometer.
3. Lift the Hand Unit from the cradle.
4. Remove the red protective dust cap from the Puff Tube.
5. When the Hand Unit is removed from the holster the two green LEDs on the front illuminate, and the pump starts, the Pulsair Intellipuff will perform a system check, when complete the display will read 'OK', refer to Section 9 for a full list of display data.

- Before using the Pulsair IntelliPuff press the Demo button to dispel any minute particles of dust or moisture which may have settled whilst the Pulsair IntelliPuff was not in use.

## 8.2 PREPARING THE PATIENT

Before using the Pulsair IntelliPuff Tonometer you should make your patient feels at ease and ensure they are in an optimum reading location, preferably with their head supported. This is because apprehension and nervousness may adversely affect the readings obtained. Follow the points outlined below to achieve this:

- Ask the patient to remove their contact lenses or spectacles if worn and to blink and breathe normally.
- Ensure that the patient is comfortable and in a relaxed position.
- To reassure the patient, you can demonstrate the procedure, using the Clear / Demo button, on the back of the patient's hand prior to taking a reading.


### Before taking a reading, you should:

- Ask the patient to blink to ensure a good and reflective tear film.
- Ensure the patient and tonometer optics are not positioned under direct lighting (i.e. spot lights or sunlight).
- Ensure the patient's eyes are fully opened. This helps to prevent squeezing, where the patient unconsciously tenses their eyelids and increases IOP.
- Throughout the reading process, you should allow the patient to blink at intervals to maintain the corneal tear film.

## 8.3 TAKING THE READING

Once the Pulsair IntelliPuff and the patient are prepared, you are ready to take a reading.

- The Pulsair IntelliPuff is set to automatically select the right eye as the first eye to be measured. If you wish to select the left eye, press the OD / OS button on the hand unit.
- Lift the Hand Unit, the pump starts and the two green LEDs illuminate.
- From a distance of about 30 cm (12 inches), look through the eyepiece and locate the patient's eye.
- Slowly move closer to the patient, maintaining alignment. Support the Pulsair IntelliPuff against your free hand and/or use the 'pop out' Forehead Rest.

- Continue to move in slowly towards the patient, two green dots will appear.
- Continue to move closer, a red reflex appears.


- Move closer. At a distance of approximately 15mm, a black cross on red or 'bow tie' image appears. Centre this image (on the central bar) and the Pulsair IntelliPuff fires automatically.

- Once you have taken a reading, remain in the operating position; wait a few seconds for the air chamber to refill. When the 'bow tie' image appears the Pulsair IntelliPuff takes a subsequent reading. When two successive readings within 1mmHg of each other are recorded, a sound may be emitted (if sounds are enabled in the User Menu Options). If successive readings of within 1mmHg of each other are not obtained, Keeler recommend taking up to four readings.
- When two consecutive readings are within 1mmHg an audible sound will be heard indicating that sufficient readings may have been taken.
- If a reading is recorded as a non-event or bad event, a long high pitch tone will be heard.
- The first reading will be the measured value; successive readings will display the running average IOP. Outlying or spurious readings will automatically be excluded from the calculation.
- At any time pressing the Review button will allow you to view the individual readings.
- If the unit does not fire, repeat step 3-7.

To measure the other eye, press the OD / OS button on the Hand Unit and repeat processes 3-7.

**Note:** For quick alignment reference, please refer to the short form instructions located at the rear of the instrument and accessed by using the pull out tab on the back left hand side.

## 9. DISPLAY EXAMPLES

STBY

**Standby**

The tonometer will display STBY when power is on.

WAIT

**System Initialisation**

The unit will display WAIT for one second while the system initialises.

OK

**OK**

When no fault is found, OK is displayed and the tonometer defaults to measure the right eye OD.

14

**Shows first reading of 14mmHg.**

The LED OD / OS indicator shows which eye the reading relates to.

14.7

**Reading to 0.1 significant figure**

If 0.1 significant figure is selected using the menu options. In this case there is no indication on the display that the reading is the first or average of more than 1.

&gt;25

**IOP greater than 25mmHg**

When a pressure (IOP) of greater than 25mmHg is detected, the unit will display >25, the puff intensity will automatically be increased from the soft puff to the normal puff level for subsequent measurements.

RUN TEST

**Self-test**

The unit will perform regular Self-tests – if a possible discrepancy in operating parameters is suspected the 'RUN TEST' message will be displayed for up to 15 seconds. To continue to use the tonometer press the OD / OS button to clear the message. The results displayed thereafter may be suspect. Refer to the User Menu Options section in this manual for guidance on running the Self-test.

ERR

**Error**

Display shows error. (Signified by a long high pitch sound).

## 10. PRINTING

The results can be printed by pressing the Print button on the Hand Unit, or if the user menu is set to automatically print by replacing the Hand Unit in the cradle.


### 10.1 SAMPLE PRINT

The reading in brackets (20) indicates a discarded value (not taken into account in the average calculations).

The Name, Date and Time fields are to be manually written in by the operator.

The average IOP is printed to one decimal place 'xx.x'.

The last four individual readings are printed to zero decimal places 'xx'.

## 11. REPLACING THE PRINTER PAPER

1. Access to the printer paper is via the Printer Cover, pull the lip on the top of the Cover and gently pull towards you to open the Printer Cover. Remove the empty paper roll.
2. Place the new roll of paper into the paper holder, making sure the free end is loose at the top of the roll, otherwise it won't print.
3. Feed the free end of paper through the gap in the Cover.
4. Close the Cover.

## 12. USER MENU OPTIONS

1. With the tonometer switched on and the Hand Unit removed, press and hold Print / Menu button for more than 3 second to enter the User Menu Options.
2. The Display will show the first User Menu Option and the current selection i.e.. PRNT ON or PRNT OFF.

- To change the User Option, press the OD OS / 'Menu Change' button once, 'toggling' the OD OS / 'Menu Change' button will cycle through the option(s).
- Pressing the Print / Menu button will move you forward to the next User Option, in this case the Buzzer Control.
- Use the OD OS / 'Menu Change' button to make your preferred selection.

**Note:** to run the Self-test, press the DEMO button, not the OD OS / 'Menu Change' button.

- Continue to repeat steps 4 and 5 until 'OK' is displayed. Your Pulsair IntelliPuff Tonometer is now ready to use with your preferred settings.

Menu Option	Display	Change Options
Printer Control	PRNT	OFF / ON
Buzzer Control	BUZ	ON / OFF
IOP Format		XX / XX.X
Full Self-test	RUN TEST	TEST / WAIT

The last choice above 'Run Test' (selected by pressing DEMO button) would start a Self-test program (around 45 seconds), the result of which should be printed.

## 13. CALIBRATION, MAINTENANCE AND INSPECTION

**Keeler recommends this routine maintenance be carried out by the user frequently to ensure safe and accurate measurement. In the event of the device being outside of the calibration tolerances, it is important to send the device back to Keeler Ltd. or your local dealer for repair and re-calibration.**

### 13.1 REGULAR INSPECTION

Inspect your power supply unit and cable for damage regularly.

Before inspecting, disconnect the power supply from the Pulsair IntelliPuff Tonometer and the mains.

If the outer insulation of the cable appears to be damaged discontinue use immediately. Contact your local dealer for a replacement.

### 13.2 GENERAL

Keep the tonometer free from dust.

If the Pulsair IntelliPuff Tonometer is to remain unused for any length of time, press the On / Off Push button switch to 'Off' and remove the power supply. Use the dust cover to protect the tonometer.

## 14. SERVICING AND CALIBRATION

Keeler recommends an annual calibration for the Tonometer. Do not modify this equipment without authorisation of the manufacturer.

This must be performed by an authorised Pulsair service centre or distributor. The unit performs a self function check when switched on and will indicate if a fault is found.

There are no user serviceable parts in this instrument. Service manuals will be available to authorised Keeler service centres and Keeler trained service personnel.

## 15. WARRANTY

Your Keeler product is guaranteed for 2 years and will be replaced, or repaired free of charge subject to the following:

- Any fault due to faulty manufacture.
- The instrument and accessories have been used in compliance with these instructions.
- Proof of purchase accompanies any claim.

**The manufacturer declines any and all responsibility and warranty coverage should the instrument be tampered with in any manner or should routine maintenance be omitted or performed in manners not in accordance with these manufacturer's instructions.**

**There are no user serviceable parts in this instrument. Any servicing or repairs should only be carried out by Keeler Ltd. or by suitably trained and authorised distributors. Service manuals will be available to authorised Keeler service centres and Keeler trained service personnel.**

## 16. SPECIFICATIONS AND ELECTRICAL RATINGS

The Keeler Pulsair IntelliPuff Tonometer is a medical electrical instrument. The instrument requires special care concerning electromagnetic compatibility (EMC). This Section describes its suitability in terms of electromagnetic compatibility of this instrument. When installing or using this instrument, please read carefully and observe what is described here.

Portable or mobile-type radio frequency communication units may have an adverse effect on this instruments, resulting in malfunctioning.

## 16.1 ELECTROMAGNETIC EMISSIONS

**Guidance and manufacturer's declaration – electromagnetic emissions**

The Keeler Pulsair IntelliPuff Tonometer is intended for use in the electromagnetic environment specified below. The customer or user should assure that it is used in such an environment.

<b>Emissions test</b>	<b>Compliance</b>	<b>Electromagnetic environment – guidance</b>
RF emissions CISPR 11	Group 1	The Keeler Pulsair IntelliPuff Tonometer uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	The Keeler Pulsair IntelliPuff Tonometer is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2		
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

## 16.2 ELECTROMAGNETIC IMMUNITY

**Guidance and manufacturer's declaration – electromagnetic immunity**

The Keeler Pulsair IntelliPuff Tonometer is intended for use in the electromagnetic environment specified below. The customer or user should assure that it is used in such an environment.

<b>Immunity test</b>	<b>IEC 55015 Test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment – guidance</b>
Electrostatic discharge (ESD). IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile.  If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst. IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge. IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.

Immunity test	IEC 55015 Test level	Compliance level	Electromagnetic environment – guidance
Voltage dips, short interruptions and voltage variations on power supply input lines. IEC 61000-4-11	$U_T = 0\%$ 0.5 cycle (0, 45, 90, 135, 180, 225, 270, 315°) $U_T = 0\%$ ; 1 cycle $U_T = 70\%$ ; 25/30 cycles (@ 0°) $U_T = 0\%$ ; 250/300 cycle	$U_T = 0\%$ 0.5 cycle (0, 45, 90, 135, 180, 225, 270, 315°) $U_T = 0\%$ ; 1 cycle $U_T = 70\%$ ; 25/30 cycles (@ 0°) $U_T = 0\%$ ; 250/300 cycle	Mains power quality should be that of a typical commercial or hospital environment.  If the user of the Keeler Pulsair Intellipuff Tonometer requires continued operation during power mains interruptions, it is recommended that the charger be powered from an uninterruptible power supply.
Power frequency (50/60 Hz) Magnetic field. IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at a level characteristic of a typical location in a typical professional healthcare facility environment.

Note:  $U_T$  is the a. c. mains voltage prior to application of the test level.

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment – guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the Keeler Pulsair Intellipuff Tonometer, including cables, than the recommended separation distances calculated from the equation applicable to the frequency of the transmitter.
<b>Recommended separation distance</b>			
Conducted RF IEC 61000-4-6	6 Vrms	6 V	$d = 1.2 \sqrt{p}$
Radiated RF IEC 61000-4-3	10 V/m 80MHz to 2.7GHz	10 V/m	$d = 1.2 \sqrt{p}$ 80MHz to 800 MHz $d = 2.3 \sqrt{p}$ 800MHz to 2.7GHz
			Where $p$ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and $d$ is the recommended separation distance in metres (m).  Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>1</sup> , should be less than the compliance level in each frequency range. <sup>2</sup>
			 Interference may occur in the vicinity of equipment marked with this symbol.

Note 1: At 80MHz and 800MHz, the higher frequency range applies.

Note 2: These guide lines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people

<sup>1</sup> Field strengths from fixed transmitters, such as base stations (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Keeler Pulsair IntelliPuff Tonometer is used exceeds the applicable RF compliance level above, the Keeler Pulsair IntelliPuff Tonometer should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orientating or relocating the Keeler Pulsair IntelliPuff Tonometer.

<sup>2</sup> Over the frequency range 150kHz to 80 MHz, field strengths should be less than 10 V/m.

### 16.3 RECOMMENDED SAFE DISTANCES

#### Recommended separation distances between portable and mobile RF communications equipment and the Keeler IntelliPuff.

The Keeler Pulsair IntelliPuff Tonometer is intended for the use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Keeler Pulsair IntelliPuff Tonometer can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Keeler Pulsair IntelliPuff Tonometer as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 230MHz $d = 1.2\sqrt{p}$	80MHz to 800MHz $d = 1.2\sqrt{p}$	800MHz to 2.7GHz $d = 2.3\sqrt{p}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where  $p$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note: 1 At 80MHz and 800MHz, the higher frequency range applies.




Note 2: These guide lines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.




## 17. TECHNICAL SPECIFICATIONS

<b>Console dimensions</b>	260 x 215 x 220mm (H x D x W)
<b>Hand Unit dimensions</b>	315 x 150 x 46mm (H x D x W)
<b>Console weight</b>	2.465Kg
<b>Hand Unit weight</b>	0.890Kg
<b>Calibrated range</b>	5mmHg to 50mmHg
<b>Repeatability (Average coefficient of variation)</b>	<5%
<b>Accuracy</b>	+/-5mmHg (95% confidence level)*
<b>Working distance</b>	20mm from surface of patient's cornea to front surface of first lens. This equates to a nominal distance of 15mm from the front of the puff tube shroud to the front surface of the patient's cornea
<b>Displayed scale</b>	4 character dot matrix scrolling
<b>Illumination system</b>	LED infra-red
<b>Length of umbilical cord</b>	2m
<b>Complies with</b>	Electrical Safety (Medical) IEC 60601-1 Electromagnetic Compatibility IEC 60601-1-2, BS EN ISO 15004-1, BS EN ISO 15004-2
<b>Power Supply Unit</b>	Switch mode, (110-240V)+/- 10% multi plug type compliant to EN 60601-1, EN 61000-6-2, EN 61000-6-3
<b>Power supply output</b>	30 VA (12V DC 2.5A)
<b>Frequency</b>	50/60 Hz

\*In a clinical study, the Pulsair IntelliPuff Tonometer appeared to slightly underestimate IOP relative to the Goldmann Tonometer at pressures above 30mmHg but these differences were not clinically significant.

## Environmental Conditions:

USE	
	
	
Shock (without packing)	10 g, duration 6 ms

STORAGE CONDITIONS		
 -10°C      55°C	 10%      95%	 700 hPa      1060 hPa
TRANSPORT CONDITIONS		
 -40°C      70°C	 10%      95%	 500 hPa      1060 hPa
Vibration, sinusoidal	10 Hz to 500 Hz: 0.5g	
Shock	30 g, duration 6 ms	
Bump	10 g, duration 6 ms	

## 18. ACCESSORIES AND SPARES

Item	Part Number
Printer paper roll	2208-L-7008
Intellipuff wall mounting kit	2414-P-7011
Pulsair Intellipuff dust cover	EP39-70304
Tonometer face shield	2415-P-7038

## 19. PACKAGING AND DISPOSAL INFORMATION

### Disposal of old electrical and electronic equipment

This symbol on the product or on its packaging and instructions indicates that this product shall not be treated as household waste.

To reduce the environmental impact of WEEE (Waste Electrical Electronic Equipment) and minimise the volume of WEEE entering landfills we encourage at product end of life that this equipment is recycled and reused.

**If you need more information on the collection reuse and recycling then please contact B2B Compliance on 01691 676124 (+44 1691 676124). (UK only).**

**Any serious incident that has occurred in relation to the device must be reported to the manufacturer and the competent authority of your Member State.**

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Copyright © Keeler Limited 2023. Publicado no Reino Unido em 2023.
Medical AG, Dornierstrasse 11
SGS United Kingdom Ltd